



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:36
In Date: 10/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:37
Out Date: 10/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105714

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5382-01	MH0GM9	1	1.13	9.75	8.39	84.2	
K5382-02	MH0GN0	2	1.14	9.96	9.19	91.3	
K5382-03	MH0GN1	3	1.12	9.97	8.44	82.7	
K5382-04	MH0GN2	4	1.12	9.95	9.58	95.8	
K5382-05	MH0GN3	5	1.16	9.52	8.44	87.1	
K5382-06	MH0GN4	6	1.13	9.76	8.73	88.1	
K5382-07	MH0GN5	7	1.14	9.95	9.1	90.4	
K5382-08	MH0GN6	8	1.13	9.97	9.2	91.3	
K5382-09	MH0GN7	9	1.16	9.96	9.54	95.2	
K5382-10	MH0GN8	10	1.14	9.88	7.76	75.7	
K5382-11	MH0GN9	11	1.14	9.84	8.68	86.7	
K5382-12	MH0GO0	12	1.12	9.76	7.86	78.0	
K5382-13	MH0GO1	13	1.14	9.92	9.39	94.0	
K5382-14	MH0GO2	14	1.12	9.65	8.36	84.9	
K5382-15	MH0GO3	15	1.12	9.95	9.28	92.4	
K5382-16	MH0GO4	16	1.14	9.7	8.62	87.4	
K5382-17	MH0GO5	17	1.14	9.87	8.66	86.1	
K5382-18	MH0GO6	18	1.14	9.98	8.92	88.0	
K5382-19	MH0GO7	19	1.15	9.62	8.82	90.6	
K5382-20	MH0GO8	20	1.16	9.84	9.02	90.6	
K5382-21	MH0GO8D	21	1.16	9.84	9.02	90.6	
K5382-22	MH0GO8S	22	1.16	9.84	9.02	90.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$