



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:47
In Date: 10/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105718

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5152-06	SU-01-101519	1	1.14	9.72	9.16	93.5	
K5152-07	SU-01-101519-1	2	1.12	9.56	9.08	94.3	
K5152-08	SU-01-101519-2	3	1.14	9.87	8.75	87.2	
K5362-01	391010741.0	4	1.16	9.92	4.09	33.4	SLUDGE SAMPLE
K5366-01	391008025.0	5	1.13	9.98	1.22	1.0	SLUDGE SAMPLE
K5367-01	391010714.0	6	1.12	9.82	1.27	1.7	SLUDGE SAMPLE
K5368-01	391010710.0	7	1.14	9.98	1.79	7.4	SLUDGE SAMPLE
K5369-01	391009706.2	8	1.14	9.81	1.46	3.7	SLUDGE SAMPLE
K5370-01	391008043.0	9	1.16	9.66	1.56	4.7	SLUDGE SAMPLE
K5371-01	391011028.0	10	1.15	9.95	1.52	4.2	SLUDGE SAMPLE
K5374-01	STOCK-PILE-1	11	1.15	9.88	8.88	88.5	
K5374-02	STOCK-PILE-1	12	1.15	9.88	8.88	88.5	
K5377-01	205-166	13	1.00	2.00	2.00	100.0	100 % SOLIDS
K5392-01	TP-22A-S1	14	1.16	9.7	8.87	90.3	
K5392-02	TP-22A-S2	15	1.15	9.67	8.92	91.2	
K5392-03	TP-22A-S	16	1.14	9.85	9.04	90.7	
K5392-05	TP-22A-D1	17	1.13	9.93	9.18	91.5	
K5392-06	TP-22A-D2	18	1.14	9.8	9.00	90.8	
K5392-07	TP-22A-D	19	1.14	9.98	9.25	91.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$