



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/17/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:48
In Date: 10/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105747

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5153-01	TR-04-101619	1	1.13	9.86	9.2	92.4	
K5153-02	TR-04-101619-1	2	1.13	9.56	8.67	89.4	
K5153-03	TR-04-101619-2	3	1.13	9.97	9.54	95.1	
K5393-01	T-6-D1-D3- (0-46)	4	1.12	9.93	8.73	86.4	
K5393-08	T-9-U1-U3- (0-28)	5	1.12	9.74	8.06	80.5	
K5393-14	T-13-V1- (0-24)	6	1.13	9.89	9.42	94.6	
K5393-17	T-16-G1- (0-27)	7	1.12	9.7	8.97	91.5	
K5393-21	T-6-D2- (46)	8	1.13	9.98	8.99	88.8	
K5393-23	T-9-U3- (28)	9	1.13	9.75	8.95	90.7	
K5393-25	T-13-V1- (24)	10	1.13	9.5	9.09	95.1	
K5393-27	T-16-G1- (27)	11	1.14	9.97	9.18	91.1	
K5394-01	T-17-F1- (0-30)	12	1.12	9.55	8.36	85.9	
K5394-04	T-15-E1- (0-20)	13	1.12	9.57	8.55	87.9	
K5394-08	T-17-F1- (30)	14	1.12	9.72	8.28	83.3	
K5394-10	T-15-E1- (20)	15	1.12	9.64	8.61	87.9	
K5400-01	55-01	16	1.12	9.61	8.99	92.7	
K5400-02	55-02	17	1.13	9.63	9.16	94.5	
K5400-03	55-03	18	1.14	9.72	9.29	95.0	
K5400-04	41-01	19	1.13	9.68	8.3	83.9	
K5400-05	41-02	20	1.13	9.58	8.33	85.2	
K5400-06	41-03	21	1.13	9.96	8.5	83.5	
K5400-07	41-04	22	1.14	9.56	8.18	83.6	
K5400-08	41-05	23	1.14	9.65	8.26	83.7	
K5400-09	41-06	24	1.15	9.73	8.33	83.7	
K5400-10	40-01	25	1.12	9.52	9.19	96.1	
K5400-11	40-02	26	1.12	9.84	9.45	95.5	
K5400-12	40-03	27	1.12	9.9	9.37	94.0	
K5400-13	40-04	28	1.13	9.72	9.42	96.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/17/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:48
In Date: 10/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105747

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5400-14	40-05	29	1.12	9.89	9.61	96.8	
K5400-15	40-06	30	1.1	9.63	9.33	96.5	
K5401-01	RBR200031-RT4221	31	1.13	9.87	8.33	82.4	
K5401-02	RBR200036-RT1488	32	1.12	9.73	8.64	87.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$