



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/21/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:52
In Date: 10/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 10/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105801

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5388-01	MH0GZ0	1	1.14	9.81	9.48	96.2	
K5388-02	MH0GZ1	2	1.14	9.5	8.8	91.6	
K5388-03	MH0GZ2	3	1.15	9.94	9.03	89.6	
K5388-04	MH0GZ3	4	1.15	9.93	9.16	91.2	
K5388-05	MH0GZ4	5	1.15	9.58	8.77	90.4	
K5388-06	MH0GZ5	6	1.14	9.91	9.06	90.3	
K5388-07	MH0GZ6	7	1.14	9.74	9.26	94.4	
K5388-08	MH0GZ7	8	1.13	9.58	8.97	92.8	
K5388-09	MH0GZ8	9	1.15	9.6	9.06	93.6	
K5388-10	MH0GZ9	10	1.16	9.57	9.24	96.1	
K5388-11	MH0HA0	11	1.16	9.74	9.12	92.8	
K5388-12	MH0HA1	12	1.16	9.6	8.98	92.7	
K5388-13	MH0HA2	13	1.14	9.55	9.00	93.5	
K5388-14	MH0HA3	14	1.14	9.7	9.32	95.6	
K5388-15	MH0HA4	15	1.14	9.52	8.65	89.6	
K5388-16	MH0HA5	16	1.15	9.61	9.09	93.9	
K5388-17	MH0HA5D	17	1.15	9.61	9.09	93.9	
K5388-18	MH0HA5S	18	1.15	9.61	9.09	93.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$