



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:07
In Date: 10/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105859

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5433-01	JLNF2	1	1.15	9.94	8.41	82.6	
K5433-02	JLNF3	2	1.12	9.82	8.39	83.6	
K5433-03	JLNF0	3	1.13	9.86	8.24	81.4	
K5433-04	JLNF0MS	4	1.13	9.86	8.24	81.4	
K5433-05	JLNF0MSD	5	1.13	9.86	8.24	81.4	
K5433-06	JLNF1	6	1.12	9.81	8.27	82.3	
K5433-07	JLNF4	7	1.12	9.63	7.86	79.2	
K5433-08	JLNF5	8	1.12	9.87	8.28	81.8	
K5433-09	JLNF6	9	1.12	9.81	6.96	67.2	
K5433-10	JLNF7	10	1.15	9.53	7.9	80.5	
K5433-11	JLNF8	11	1.13	9.89	7.71	75.1	
K5433-12	JLNG1	12	1.12	9.6	7.56	75.9	
K5433-13	JLNF9	13	1.13	9.76	6.83	66.0	
K5433-14	JLNG0	14	1.14	9.61	7.54	75.6	
K5433-15	JLNG2	15	1.14	9.8	8.62	86.4	
K5433-16	JLNG3	16	1.13	9.9	8.29	81.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$