



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:31
In Date: 10/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:43
Out Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105860

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5434-01	JLNG4	1	1.14	9.89	7.88	77.0	
K5434-02	JLNG4MS	2	1.14	9.89	7.88	77.0	
K5434-03	JLNG4MSD	3	1.14	9.89	7.88	77.0	
K5434-04	JLNG5	4	1.12	9.73	6.61	63.8	
K5434-05	JLNG6	5	1.12	9.83	6.3	59.5	
K5434-06	JLNG7	6	1.13	9.98	5.6	50.5	
K5434-07	JLNG8	7	1.12	9.61	7.67	77.1	
K5434-08	JLNG9	8	1.12	9.86	3.51	27.3	
K5434-09	JLNH0	9	1.12	9.89	4.54	39.0	
K5434-10	JLNH1	10	1.14	9.72	4.35	37.4	
K5434-11	JLNH2	11	1.13	9.65	7.88	79.2	
K5434-12	JLNH3	12	1.12	9.86	4.35	37.0	
K5434-13	JLNH4	13	1.14	9.97	7.34	70.2	
K5434-14	JLNH5	14	1.14	9.83	7.08	68.4	
K5434-15	JLNH6	15	1.13	9.61	7.29	72.6	
K5434-16	JLNH7	16	1.13	9.8	7.77	76.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$