



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/24/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:19
In Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105879

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5145-04	HD-02-102319	135	1.14	9.79	8.41	84.0	
K5145-05	HD-02-102319	136	1.13	9.78	8.53	85.5	
K5148-09	OR-02-102319	137	1.15	9.65	8.78	89.8	
K5148-10	OR-02-102319-1	138	1.14	9.66	8.75	89.3	
K5148-11	OR-02-102319-2	139	1.13	9.7	8.5	86.0	
K5489-03	MAIN-STOCKPILE-SLP	1	1.12	9.81	9.14	92.3	
K5497-01	170-48-(0-1)	2	1.13	9.61	8.57	87.7	
K5497-02	170-48-(1)	3	1.16	9.89	8.77	87.2	
K5497-03	169-90-(0-1.2)	4	1.16	9.66	8.79	89.8	
K5497-04	169-90-(1.2)	5	1.12	9.9	8.78	87.2	
K5497-05	169-32-(0-1)	6	1.13	9.74	8.82	89.3	
K5497-06	169-32-(1)	7	1.15	9.55	8.58	88.5	
K5497-07	168-74-(0-1)	8	1.12	9.65	8.16	82.5	
K5497-08	168-74-(1)	9	1.12	9.51	8.03	82.4	
K5497-09	168-16-(0-1)	10	1.14	9.67	6.06	57.7	
K5497-10	168-16-(1)	11	1.13	9.77	8.23	82.2	
K5497-11	167-58-(0-1)	12	1.15	9.9	8.96	89.3	
K5497-12	167-58-(1)	13	1.12	9.72	8.6	87.0	
K5497-13	167-00-(0-1)	14	1.14	9.5	6.4	62.9	
K5497-14	167-00-(1)	15	1.12	9.8	6.72	64.5	
K5497-15	166-42-(0-1)	16	1.12	9.73	8.24	82.7	
K5497-16	166-42-(1)	17	1.13	9.81	8.43	84.1	
K5497-17	188-46-(0-1)	18	1.15	9.76	9.18	93.3	
K5497-18	188-46-(1)	19	1.14	9.98	9.23	91.5	
K5497-19	189-04-(0-1.2)	20	1.14	9.96	9.26	92.1	
K5497-20	189-04-(1.2)	21	1.15	9.53	8.88	92.2	
K5498-01	189-62-(0-1)	22	1.13	9.86	9.2	92.4	
K5498-02	189-62-(1)	23	1.14	9.98	9.45	94.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/24/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:19
In Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105879

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5498-03	190-20-(0-1)	24	1.14	9.8	8.79	88.3	
K5498-04	190-20-(1)	25	1.14	9.6	8.58	87.9	
K5498-05	190-78-(0-1)	26	1.14	9.55	8.29	85.0	
K5498-06	190-78-(1)	27	1.14	9.92	8.76	86.8	
K5498-07	196-00-(0-1)	28	1.14	9.81	8.77	88.0	
K5498-08	196-00-(1)	29	1.15	9.5	8.44	87.3	
K5498-09	196-58-(0-1)	30	1.16	9.55	7.57	76.4	
K5498-10	196-58-(1)	31	1.13	9.85	7.63	74.5	
K5498-11	191-36-(0-1)	32	1.12	9.61	8.00	81.0	
K5498-12	191-36-(1)	33	1.14	9.73	7.92	78.9	
K5498-13	191-94-(0-1)	34	1.15	9.68	8.26	83.4	
K5498-14	194-94-(1)	35	1.14	9.54	8.26	84.8	
K5498-15	191-52-(0-1)	36	1.13	9.82	7.87	77.6	
K5498-16	191-52-(1)	37	1.13	9.96	8.00	77.8	
K5498-17	193-10-(0-1)	38	1.12	9.9	7.54	73.1	
K5498-18	193-10-(1)	39	1.12	9.68	7.13	70.2	
K5498-19	193-68-(0-1)	40	1.15	9.82	7.35	71.5	
K5498-20	193-68-(1)	41	1.12	9.58	7.17	71.5	
K5499-01	193-26-(0-1)	42	1.13	9.98	7.53	72.3	
K5499-02	193-26-(1)	43	1.12	9.54	6.8	67.5	
K5499-03	194-84-(0-1)	44	1.14	9.92	7.95	77.6	
K5499-04	194-84-(1)	45	1.12	9.66	7.69	76.9	
K5499-05	195-42-(0-1)	46	1.13	9.68	6.57	63.6	
K5499-06	195-42-(1)	47	1.12	9.76	6.58	63.2	
K5499-07	197-16-(0-1)	48	1.15	9.55	6.32	61.5	
K5499-08	197-16-(1)	49	1.13	9.5	5.89	56.9	
K5499-09	197-74-(0-1)	50	1.14	9.58	7.97	80.9	
K5499-10	197-74-(1)	51	1.14	9.73	8.48	85.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/24/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:19
In Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105879

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5499-11	198-27-(0-1.8)	52	1.14	9.64	8.58	87.5	
K5499-12	198-27-(1.8)	53	1.13	9.5	8.47	87.7	
K5499-13	707-58-(0-1)	54	1.13	9.78	9.27	94.1	
K5499-14	707-58-(1)	55	1.13	9.78	9.21	93.4	
K5499-15	705-81-(0-1)	56	1.13	9.6	6.87	67.8	
K5499-16	705-81-(1)	57	1.12	9.98	7.08	67.3	
K5499-17	704-04-(0-1)	58	1.16	9.7	6.43	61.7	
K5499-18	704-04-(1)	59	1.15	9.82	6.87	66.0	
K5499-19	702-27-(0-1)	60	1.15	9.57	8.63	88.8	
K5499-20	702-27-(1)	61	1.14	9.62	8.74	89.6	
K5500-01	802-50-(0-1.1)	62	1.16	9.76	8.98	90.9	
K5500-02	802-50-(1.1)	63	1.14	9.78	8.88	89.6	
K5500-03	802-00-(0-2)	64	1.15	9.63	8.72	89.3	
K5500-04	802-00-(2)	65	1.14	9.93	8.88	88.1	
K5500-05	801-50-(0-1.1)	66	1.16	9.62	8.89	91.4	
K5500-06	801-50-(1.1)	67	1.14	9.8	8.96	90.3	
K5500-07	804-00-(0-1.8)	68	1.14	9.51	8.44	87.2	
K5500-08	804-00-(1.8)	69	1.14	9.94	8.83	87.4	
K5500-09	800-50-(0-2)	70	1.13	9.78	8.67	87.2	
K5500-10	800-50-(2)	71	1.13	9.87	8.65	86.0	
K5501-01	187-88-(0-2)	72	1.15	9.53	8.69	90.0	
K5501-02	187-88-(2)	73	1.14	9.97	9.18	91.1	
K5501-03	187-30-(0-2)	74	1.15	9.89	9.03	90.2	
K5501-04	187-30-(2)	75	1.14	9.98	9.01	89.0	
K5501-05	186-72-(0-2.2)	76	1.13	9.76	9.28	94.4	
K5501-06	186-72-(2.2)	77	1.12	9.7	9.18	93.9	
K5501-07	186-14-(0-2.1)	78	1.14	9.5	9.03	94.4	
K5501-08	186-14-(2.1)	79	1.13	9.62	9.18	94.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/24/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:19
In Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105879

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5501-09	185-56- (0-2.2)	80	1.12	9.72	9.22	94.2	
K5501-10	185-56- (2.2)	81	1.12	9.6	9.05	93.5	
K5501-11	184-98- (0-2.2)	82	1.12	9.57	9.04	93.7	
K5501-12	184-98- (2.2)	83	1.13	9.64	9.11	93.8	
K5501-13	184-40- (0-2)	84	1.14	9.91	9.46	94.9	
K5501-14	184-40- (2)	85	1.12	9.64	9.18	94.6	
K5501-15	183-82- (0-2.1)	86	1.13	9.57	8.98	93.0	
K5501-16	183-82- (2.1)	87	1.14	9.97	9.41	93.7	
K5501-17	183-24- (0-2.1)	88	1.13	9.64	9.16	94.4	
K5501-18	183-24- (2.1)	89	1.14	9.62	9.14	94.3	
K5501-19	182-66- (0-2)	90	1.14	9.85	8.65	86.2	
K5501-20	182-66- (2)	91	1.14	9.82	8.59	85.8	
K5502-01	182-08- (0-2)	92	1.14	9.92	8.9	88.4	
K5502-02	182-08- (2)	93	1.14	9.67	8.43	85.5	
K5502-03	181-50- (0-2)	94	1.16	9.91	8.06	78.9	
K5502-04	181-50- (2)	95	1.14	9.67	8.17	82.4	
K5502-05	180-92- (0-2)	96	1.14	9.53	8.55	88.3	
K5502-06	180-92- (2)	97	1.15	9.61	8.67	88.9	
K5502-07	180-34- (0-2.1)	98	1.14	9.89	8.55	84.7	
K5502-08	180-34- (2.1)	99	1.16	9.62	8.3	84.4	
K5502-09	179-76- (0-2.1)	100	1.15	9.6	7.82	78.9	
K5502-10	179-76- (2.1)	101	1.16	9.87	8.02	78.8	
K5502-11	179-18- (0-2.1)	102	1.14	9.58	8.39	85.9	
K5502-12	179-18- (2.1)	103	1.16	9.89	8.72	86.6	
K5502-13	178-60- (0-2.1)	104	1.13	9.63	8.6	87.9	
K5502-14	178-60- (2.1)	105	1.14	9.91	8.73	86.5	
K5502-15	178-02- (0-2.3)	106	1.15	9.57	8.27	84.6	
K5502-16	178-02- (2.3)	107	1.14	9.55	8.18	83.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/24/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:19
In Date: 10/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105879

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5502-17	177-44- (0-2.5)	108	1.15	9.59	8.28	84.5	
K5502-18	177-44- (2.5)	109	1.13	9.58	8.19	83.6	
K5502-19	176-86- (0-2.5)	110	1.15	9.72	8.35	84.0	
K5502-20	176-86- (2.5)	111	1.14	9.98	8.61	84.5	
K5503-01	176-28- (0-2.5)	112	1.16	9.9	8.67	85.9	
K5503-02	176-28- (2.5)	113	1.12	9.85	8.49	84.4	
K5503-03	175-70- (0-2.5)	114	1.14	9.5	8.25	85.0	
K5503-04	175-70- (2.5)	115	1.13	9.56	8.21	84.0	
K5503-05	175-12- (0-2)	116	1.15	9.73	8.72	88.2	
K5503-06	175-12- (2)	117	1.15	9.62	8.45	86.2	
K5503-07	174-54- (0-2)	118	1.14	9.89	8.83	87.9	
K5503-08	174-54- (2)	119	1.16	9.5	8.46	87.5	
K5503-09	173-96- (0-2)	120	1.14	9.59	8.25	84.1	
K5503-10	173-96- (2)	121	1.12	9.65	8.26	83.7	
K5503-11	173-38- (0-2)	122	1.15	9.53	7.32	73.6	
K5503-12	173-38- (2)	123	1.14	9.78	7.33	71.6	
K5503-13	172-80- (0-2)	124	1.13	9.57	7.84	79.5	
K5503-14	172-80- (2)	125	1.12	9.71	8.43	85.1	
K5503-15	172-22- (0-2)	126	1.13	9.55	7.49	75.5	
K5503-16	172-22- (2)	127	1.14	9.82	7.46	72.8	
K5503-17	171-64- (0-1.8)	128	1.15	9.72	8.3	83.4	
K5503-18	171-64- (1.8)	129	1.13	9.53	8.06	82.5	
K5503-19	171-06- (0-1.9)	130	1.13	9.56	8.33	85.4	
K5503-20	171-06- (1.9)	131	1.12	9.89	8.08	79.4	
K5510-01	TOPSOIL	132	1.14	9.6	8.51	87.1	
K5510-02	STOCK-PILE	133	1.17	9.98	9.12	90.2	
K5512-01	CHRT-26706	140	1.16	9.81	9.38	95.0	
K5515-01	Y-EDISON	134	1.16	9.79	8.77	88.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$