



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105914

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5153-05	TR-05-102419	1	1.16	9.64	8.68	88.7	
K5153-06	TR-05-102419	2	1.15	9.54	8.58	88.6	
K5525-01	391017139	3	1.16	9.86	1.51	4.0	SLUDGE SAMPLE
K5526-01	391018700	4	1.15	9.64	3.21	24.3	SLUDGE SAMPLE
K5533-01	P001-CHIP001-01	5	1.00	2.00	2.00	100.0	PAINTCHIPS 100 % SOLIDS
K5533-02	P001-CHIP002-01	6	1.00	2.00	2.00	100.0	PAINTCHIPS 100 % SOLIDS
K5534-01	OIL-E	7	1.00	2.00	2.00	100.0	OIL SAMPLE 100%
K5534-03	OIL-A-C	8	1.00	2.00	2.00	100.0	OIL SAMPLE 100%
K5534-05	OILY-DEBRIS	9	1.00	2.00	2.00	100.0	OILY DEBRIS 100% SOLIDS
K5539-01	160-62-(0-1)	10	1.15	9.97	8.83	87.1	
K5539-02	160-62-(1)	11	1.13	9.52	8.28	85.2	
K5539-03	161-16-(0-1)	12	1.14	9.87	8.66	86.1	
K5539-04	161-16-(1)	13	1.15	9.98	8.81	86.7	
K5539-05	161-78-(0-1.1)	14	1.14	9.73	8.28	83.1	
K5539-06	161-78-(1.1)	15	1.14	9.71	7.98	79.8	
K5539-07	162-36-(0-1)	16	1.15	9.65	8.52	86.7	
K5539-08	162-36-(1)	17	1.14	9.63	8.42	85.7	
K5539-09	162-94-(0-1)	18	1.13	9.98	8.39	82.0	
K5539-10	162-94-(1)	19	1.14	9.73	8.25	82.8	
K5539-11	163-52-(0-1)	20	1.15	9.97	8.95	88.4	
K5539-12	163-52-(1)	21	1.14	9.65	8.68	88.6	
K5539-13	164-10-(0-1)	22	1.15	9.71	8.12	81.4	
K5539-14	164-10-(1)	23	1.16	9.56	7.53	75.8	
K5539-15	164-68-(0-1)	24	1.14	9.6	7.58	76.1	
K5539-16	164-68-(1)	25	1.14	9.76	7.6	74.9	
K5539-17	165-26-(0-1)	26	1.15	9.89	8.11	79.6	
K5539-18	165-26-(1)	27	1.14	9.7	7.82	78.0	
K5539-19	165-84-(0-1)	28	1.13	9.5	6.61	65.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105914

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5539-20	165-84-(1)	29	1.13	9.97	8.37	81.9	
K5540-01	154-82-(0-2)	30	1.16	9.62	6.59	64.2	
K5540-02	154-82-(2)	31	1.15	9.6	6.51	63.4	
K5540-03	155-40-(0-1.1)	32	1.15	9.7	7.28	71.7	
K5540-04	155-40-(1.1)	33	1.14	9.67	7.34	72.7	
K5540-05	155-98-(0-1.1)	34	1.13	9.55	6.48	63.5	
K5540-06	155-98-(1.1)	35	1.13	9.7	6.88	67.1	
K5540-07	156-56-(0-1.3)	36	1.13	9.9	5.3	47.5	
K5540-08	156-56-(1.3)	37	1.15	9.75	6.5	62.2	
K5540-09	157-14-(0-1.6)	38	1.13	9.64	6.99	68.9	
K5540-10	157-14-(1.6)	39	1.15	9.84	7.21	69.7	
K5540-11	157-72-(0-1.1)	40	1.14	9.59	6.98	69.1	
K5540-12	157-72-(1.1)	41	1.15	9.93	6.97	66.3	
K5540-13	158-30-(0-1)	42	1.15	9.52	6.85	68.1	
K5540-14	158-30-(1)	43	1.14	9.98	7.64	73.5	
K5540-15	158-88-(0-1)	44	1.13	9.82	7.89	77.8	
K5540-16	158-88-(1)	45	1.13	9.98	8.27	80.7	
K5540-17	159-46-(0-1)	46	1.13	9.61	8.15	82.8	
K5540-18	159-46-(1)	47	1.13	9.56	7.91	80.4	
K5540-19	160-04-(0-1)	48	1.14	9.75	8.54	85.9	
K5540-20	160-04-(1)	49	1.14	9.96	8.8	86.8	
K5541-01	154-24-(0-2)	50	1.14	9.98	8.47	82.9	
K5541-02	154-24-(2)	51	1.13	9.6	7.99	81.0	
K5541-03	153-63-(0-2)	52	1.14	9.73	7.86	78.2	
K5541-04	153-63-(2)	53	1.15	9.98	8.38	81.9	
K5541-05	153-08-(0-1.1)	54	1.15	9.55	7.35	73.8	
K5541-06	153-08-(1.1)	55	1.13	9.79	7.5	73.6	
K5541-07	152-50-(0-1)	56	1.13	9.83	8.94	89.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105914

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5541-08	152-50-(1)	57	1.16	9.77	8.9	89.9	
K5541-09	51-00-(0-1)	58	1.14	9.54	7.56	76.4	
K5541-10	51-00-(1)	59	1.14	9.65	8.02	80.8	
K5541-11	17-20-(0-2.5)	60	1.14	9.66	8.26	83.6	
K5541-12	17-20-(2.5)	61	1.15	9.76	8.53	85.7	
K5541-13	15-70-(0-2)	62	1.12	9.64	8.01	80.9	
K5541-14	15-70-(2)	63	1.13	9.67	8.2	82.8	
K5541-15	18-12-(0-2.5)	64	1.14	9.85	8.77	87.6	
K5541-16	18-12-(2.5)	65	1.13	9.77	8.59	86.3	
K5541-17	18-48-(0-2.5)	66	1.13	9.62	8.62	88.2	
K5541-18	18-48-(2.5)	67	1.14	9.76	8.68	87.5	
K5541-19	52-77-(0-1.8)	68	1.16	9.9	9.26	92.7	
K5541-20	52-77-(1.8)	69	1.15	9.98	9.47	94.2	
K5542-01	54-54-(0-1.2)	70	1.14	9.59	8.98	92.8	
K5542-02	54-54-(1.2)	71	1.16	9.97	9.48	94.4	
K5542-03	199-48-(0-1.2)	72	1.15	9.84	8.39	83.3	
K5542-04	199-48-(1.2)	73	1.14	9.96	8.55	84.0	
K5542-05	200-06-(0-1.2)	74	1.14	9.55	8.07	82.4	
K5542-06	200-06-(1.2)	75	1.14	9.68	8.26	83.4	
K5542-07	198-90-(0-2)	76	1.14	9.59	8.27	84.4	
K5542-08	198-90-(2)	77	1.12	9.56	8.34	85.5	
K5542-09	200-64-(0-1)	78	1.11	9.88	8.02	78.8	
K5542-10	200-64-(1)	79	1.15	9.66	8.26	83.5	
K5542-11	201-22-(0-1)	80	1.14	9.54	7.86	80.0	
K5542-12	201-22-(1)	81	1.14	9.74	8.05	80.3	
K5542-13	201-80-(0-1)	82	1.15	9.89	8.5	84.1	
K5542-14	201-80-(1)	83	1.15	9.9	8.52	84.2	
K5542-15	202-38-(0-1.5)	84	1.14	9.55	8.44	86.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 10/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105914

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5542-16	202-38-(1.5)	85	1.15	9.55	8.51	87.6	
K5542-17	608-60-(0-1.5)	86	1.14	9.67	8.19	82.6	
K5542-18	608-60-(1.5)	87	1.15	9.69	8.17	82.2	
K5543-01	R-1	88	1.16	9.58	9.11	94.4	
K5543-02	R-1A	89	1.15	9.87	9.38	94.4	
K5544-01	N-1	90	1.16	9.56	9.19	95.6	
K5544-02	N-1A	91	1.14	9.52	9.19	96.1	
K5545-01	C-1	92	1.14	9.94	8.89	88.1	
K5545-02	C-1A	93	1.15	9.63	8.82	90.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$