



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/29/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:51
In Date: 10/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105969

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5149-08	EO-01-102819	18	1.14	9.68	8.8	89.7	
K5149-09	EO-01-102819-1	19	1.14	9.65	8.47	86.1	
K5149-10	EO-01-102819-2	20	1.14	9.6	8.39	85.7	
K5149-11	EO-01-102819-3	21	1.14	9.85	8.49	84.4	
K5427-20	445X-7B	1	1.13	9.95	8.98	89.0	
K5428-03	445X-8B	2	1.15	9.85	7.89	77.5	
K5428-06	445X-9B	3	1.15	9.75	7.64	75.5	
K5431-18	445X-DUP-2	4	1.13	9.63	7.57	75.8	
K5576-02	CSA-2-1-4-2	5	1.13	9.7	9.12	93.2	
K5576-04	CSA-2-1-5-2	6	1.14	9.92	9.2	91.8	
K5576-06	CSA-2-1-6-2	7	1.12	9.65	8.87	90.9	
K5576-08	CSA-2-1-7-2	8	1.14	9.83	9.36	94.6	
K5576-10	CSA-2-1-8-2	9	1.12	9.67	9.08	93.1	
K5576-14	CSA-2-1-3	10	1.13	9.71	9.03	92.1	
K5576-15	CSA-2-1-4	11	1.15	9.71	9.12	93.1	
K5576-16	CSA-2-1-5	12	1.14	9.74	8.98	91.2	
K5576-17	CSA-2-1-6	13	1.14	9.53	8.73	90.5	
K5576-18	CSA-2-1-7	14	1.14	9.63	9.22	95.2	
K5576-19	CSA-2-1-8	15	1.14	9.83	8.6	85.8	
K5584-01	445X-9A- (0-2)	16	1.14	9.62	7.67	77.0	
K5584-02	445X-9A- (2-4)	17	1.12	9.78	7.93	78.6	
K5587-01	CSA-2-1-9-2	22	1.14	9.78	8.81	88.8	
K5587-02	CSA-2-1-9-4	23	1.15	9.79	9.06	91.6	
K5587-04	CSA-2-1-10-2	24	1.15	9.69	8.93	91.1	
K5587-05	CSA-2-1-10-4	25	1.14	9.62	8.91	91.6	
K5587-07	CSA-2-1-11-2	26	1.14	9.76	9.02	91.4	
K5587-08	CSA-2-1-11-4	27	1.15	9.95	8.62	84.9	
K5587-10	CSA-2-1-12-2	28	1.14	9.56	8.38	86.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/29/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:51
In Date: 10/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105969

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5587-11	CSA-2-1-12-4	29	1.14	9.94	8.56	84.3	
K5587-13	CSA-2-1-9	30	1.14	9.56	8.67	89.4	
K5587-14	CSA-2-1-10	31	1.14	9.75	8.89	90.0	
K5587-15	CSA-2-1-11	32	1.15	9.9	8.97	89.4	
K5587-16	CSA-2-1-12	33	1.14	9.79	8.79	88.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$