



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/30/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 13:39
In Date: 10/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:39
Out Date: 10/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105985

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5596-01	GAQA2	1	1.13	9.98	8.22	80.1	
K5596-02	GAQA3	2	1.13	9.68	7.89	79.1	
K5596-03	GAQA4	3	1.12	9.61	7.61	76.4	
K5596-04	GAQA5	4	1.12	9.61	7.41	74.1	
K5596-05	GAQA6	5	1.12	9.98	7.75	74.8	
K5596-06	GAQA7	6	1.15	9.97	8.46	82.9	
K5596-07	GAQA8	7	1.15	9.84	9.18	92.4	
K5596-08	GAQA9	8	1.14	9.91	8.02	78.4	
K5596-09	GAQB0	9	1.14	9.51	7.99	81.8	
K5596-10	GAQB1	10	1.14	9.73	7.92	78.9	
K5596-11	GAQB2	11	1.15	9.72	8.06	80.6	
K5596-12	GAQB3	12	1.12	9.94	7.97	77.7	
K5596-13	GAQB4	13	1.12	9.54	7.78	79.1	
K5596-14	GAQB5	14	1.13	9.67	8.06	81.1	
K5596-15	GAQB6	15	1.14	9.66	7.88	79.1	
K5596-16	GAQB7	16	1.12	9.65	7.43	74.0	
K5596-17	GAQB8	17	1.14	9.71	7.8	77.7	
K5596-18	GAQB9	18	1.13	9.68	8.00	80.4	
K5596-19	GAQC0	19	1.14	9.57	7.78	78.8	
K5596-20	GAQC1	20	1.13	9.55	7.83	79.6	
K5596-21	GAQC1MS	21	1.13	9.55	7.83	79.6	
K5596-22	GAQC1MSD	22	1.13	9.55	7.83	79.6	
K5596-23	VHBLK01	23	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$