



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/4/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:19
In Date: 11/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106065

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5644-01	391030724	1	1.16	9.9	1.47	3.5	SLUDGE SAMPLE
K5655-05	SVOC-GPC-BLANK	2	1.00	2.00	2.00	100.0	VHBLK
K5655-06	PEST-GPC-BLANK	3	1.00	2.00	2.00	100.0	VHBLK
K5655-07	PEST-GPC-BLANK-SPIKE	4	1.00	2.00	2.00	100.0	VHBLK
K5655-08	PCB-GPC-BLANK	5	1.00	2.00	2.00	100.0	VHBLK
K5655-09	PCB-GPC-BLANK-SPIKE	6	1.00	2.00	2.00	100.0	VHBLK
K5655-10	SVOC-GPC2-BLANK	7	1.00	2.00	2.00	100.0	VHBLK
K5655-11	PEST-GPC2-BLANK	8	1.00	2.00	2.00	100.0	VHBLK
K5655-12	PEST-GPC2-BLANK-SPIKE	9	1.00	2.00	2.00	100.0	VHBLK
K5655-13	PCB-GPC2-BLANK	10	1.00	2.00	2.00	100.0	VHBLK
K5655-14	PCB-GPC2-BLANK -SPIKE	11	1.00	2.00	2.00	100.0	VHBLK
K5657-01	TP	32	1.14	9.84	6.94	66.7	
K5658-01	TP-30-S1	33	1.15	9.86	8.92	89.2	
K5658-02	TP-30-S2	34	1.14	9.67	9.38	96.6	
K5658-03	TP-30-S	35	1.16	9.79	9.4	95.5	
K5658-05	TP-30-D1	36	1.13	9.7	9.16	93.7	
K5658-06	TP-30-D2	37	1.15	9.89	8.99	89.7	
K5658-07	TP-30-D	38	1.15	9.78	9.24	93.7	
K5660-01	(A) -DT4-1- (5)	12	1.14	9.64	8.58	87.5	
K5660-02	(A) -DT4-2- (5)	13	1.15	9.62	7.92	79.9	
K5660-03	(A) -DT4-3- (5)	14	1.15	9.85	8.46	84.0	
K5660-04	(A) -DT4-4- (5)	15	1.16	9.67	8.04	80.8	
K5660-05	(A) -DT4-5- (5)	16	1.15	9.76	7.99	79.4	
K5660-06	(A) -DT4-6- (5)	17	1.15	9.86	8.37	82.9	
K5660-07	(A) -DT4-7- (5)	18	1.14	9.71	8.25	83.0	
K5660-08	(A) -DT4-8- (5)	19	1.15	9.91	8.45	83.3	
K5660-09	(A) -DT4-9- (10)	20	1.15	9.52	6.83	67.9	
K5660-10	(A) -DT4-10- (10)	21	1.14	9.93	6.9	65.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/4/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:19
In Date: 11/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106065

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5660-11	(A)-DT4-11-(10)	22	1.15	9.52	6.89	68.6	
K5660-12	(A)-DT4-12-(10)	23	1.14	9.55	7.19	71.9	
K5660-13	(A)-DT4-13-(10)	24	1.14	9.68	6.73	65.5	
K5660-14	(A)-DT4-14-(10)	25	1.16	9.78	7.62	74.9	
K5660-15	(A)-DT4-15-(10)	26	1.16	9.69	7.58	75.3	
K5660-16	(A)-DT4-16-(10)	27	1.15	9.96	7.23	69.0	
K5661-01	T-1	28	1.15	9.89	8.77	87.2	
K5661-02	T-2	29	1.15	9.93	8.79	87.0	
K5662-01	S-1	30	1.14	9.88	8.78	87.4	
K5662-02	S-1A	31	1.14	9.67	8.63	87.8	
K5676-01	SU-02-110119	39	1.15	9.69	8.2	82.6	
K5676-02	SU-02-110119	40	1.14	9.81	8.36	83.3	
K5676-03	SU-02-110119-1	41	1.16	9.97	8.85	87.3	
K5676-04	SU-02-110119-2	42	1.16	9.65	8.46	86.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$