



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 11/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 11/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106084

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5677-01	MJLM53	1	1.15	9.85	9.63	97.5	
K5677-02	MJLM57	2	1.14	9.63	9.51	98.6	
K5677-03	MJLM61	3	1.15	9.81	9.54	96.9	
K5677-04	MJLM65	4	1.14	9.9	9.54	95.9	
K5677-05	MJLM69	5	1.14	9.68	9.18	94.1	
K5677-06	MJLM73	6	1.14	9.84	9.06	91.0	
K5677-07	MJLM77	7	1.14	9.6	9.17	94.9	
K5677-08	MJLM81	8	1.12	9.94	9.54	95.5	
K5677-09	MJLM85	9	1.14	9.55	9.08	94.4	
K5677-10	MJLM89	10	1.14	9.8	9.61	97.8	
K5677-11	MJLM93	11	1.14	9.96	9.8	98.2	
K5677-12	MJLM97	12	1.14	9.65	9.51	98.4	
K5677-13	MJLM98	13	1.14	9.74	9.2	93.7	
K5677-14	MJLM99	14	1.13	9.69	9.65	99.5	
K5677-15	MJLMA0	15	1.15	9.53	9.14	95.3	
K5677-16	MJLMA1	16	1.15	9.53	9.42	98.7	
K5677-17	MJLMA1D	17	1.15	9.53	9.42	98.7	
K5677-18	MJLMA1S	18	1.15	9.53	9.42	98.7	
K5677-19	MJLMK4	19	1.15	9.83	9.7	98.5	
K5677-20	MJLMK6	20	1.14	9.65	9.52	98.5	
K5677-21	MJLMM1	21	1.14	9.87	9.86	99.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$