



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 08:10
In Date: 11/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106128

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5668-01	HD-01-110619	54	1.15	9.73	8.76	88.7	
K5668-02	HD-01-110619	55	1.19	9.66	8.44	85.6	
K5668-03	HD-01-110619-1	56	1.12	9.59	8.09	82.3	
K5675-06	TR-04-110619	57	1.15	9.96	9.12	90.5	
K5675-07	TR-04-110619	58	1.18	9.98	9.45	94.0	
K5676-06	SU-04-110619	59	1.15	9.73	8.76	88.7	
K5676-07	SU-04-110619	60	1.11	9.86	8.98	89.9	
K5676-08	SU-04-110619-1	61	1.18	9.72	8.81	89.3	
K5676-09	SU-04-110619-2	62	1.18	9.82	9.01	90.6	
K5710-01	TP- 31-S1	1	1.14	9.64	7.95	80.1	
K5710-02	TP- 31-S2	2	1.15	9.94	8.66	85.4	
K5710-03	TP- 31-S	3	1.18	9.72	8.6	86.9	
K5710-05	TP- 32-S1	4	1.18	9.88	9.03	90.2	
K5710-06	TP- 32-S2	5	1.16	9.93	8.76	86.7	
K5710-07	TP- 32-S	6	1.13	9.93	8.67	85.7	
K5710-09	TP- 32-D1	7	1.11	9.5	8.41	87.0	
K5710-10	TP- 32-D2	8	1.15	9.9	9.35	93.7	
K5710-11	TP- 32-D	9	1.12	9.59	8.61	88.4	
K5712-01	SB-2A	10	1.15	9.88	9.36	94.0	
K5712-02	SB-3A	11	1.15	9.97	9.33	92.7	
K5713-02	TT-016-IDWSO-20191105	12	1.13	9.88	7.3	70.5	
K5713-03	TT-017-IDWSO-20191105	13	1.16	9.79	7.71	75.9	
K5717-01	VC-TP	63	1.15	9.93	9.36	93.5	
K5717-02	VC-TP	64	1.1	9.67	9.15	93.9	
K5717-03	VC-TP	65	1.15	9.93	9.36	93.5	
K5722-01	4-A	14	1.12	9.68	8.61	87.5	
K5722-02	4-B	15	1.17	9.69	8.58	87.0	
K5722-03	4-C	16	1.11	9.88	8.68	86.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 08:10
In Date: 11/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106128

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5722-04	4-D	17	1.15	9.52	8.5	87.8	
K5722-05	5-A	18	1.15	9.94	9.21	91.7	
K5722-06	5-B	19	1.1	9.81	9.00	90.7	
K5722-07	5-C	20	1.14	9.72	8.93	90.8	
K5722-08	6-A	21	1.18	9.6	8.74	89.8	
K5722-09	6-B	22	1.16	9.76	8.86	89.5	
K5722-10	6-C	23	1.18	9.9	8.97	89.3	
K5722-11	1-A	24	1.15	9.68	9.62	99.3	
K5722-12	1-B	25	1.1	9.71	9.64	99.2	
K5722-13	1-C	26	1.18	9.64	9.51	98.5	
K5722-14	2-A	27	1.19	9.8	9.02	90.9	
K5722-15	2-B	28	1.16	9.92	9.14	91.1	
K5722-16	2-C	29	1.17	9.74	8.95	90.8	
K5722-17	3-A	30	1.14	9.61	8.63	88.4	
K5722-18	3-B	31	1.11	9.94	8.91	88.3	
K5722-19	3-C	32	1.18	9.76	8.81	88.9	
K5722-20	3-D	33	1.19	9.89	8.96	89.3	
K5723-01	7-A	34	1.15	9.95	8.47	83.2	
K5723-02	7-B	35	1.14	9.89	8.92	88.9	
K5723-03	7-C	36	1.11	9.8	9.00	90.8	
K5723-04	8-A	37	1.11	9.96	9.13	90.6	
K5723-05	8-B	38	1.15	9.67	8.8	89.8	
K5723-06	8-C	39	1.18	9.7	8.8	89.4	
K5723-07	8-D	40	1.11	9.97	9.03	89.4	
K5723-08	9-A	41	1.15	9.74	8.99	91.3	
K5723-09	9-B	42	1.16	9.96	9.24	91.8	
K5723-10	9-C	43	1.18	9.69	8.99	91.8	
K5723-11	9-D	44	1.12	9.85	9.05	90.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 08:10
In Date: 11/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106128

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5723-12	10-A	45	1.18	9.97	9.45	94.1	
K5723-13	10-B	46	1.18	9.72	9.15	93.3	
K5723-14	10-C	47	1.18	9.9	9.3	93.1	
K5723-15	10-D	48	1.16	9.91	9.41	94.3	
K5723-16	4-E	49	1.15	9.98	9.18	90.9	
K5725-01	SO-1	50	1.15	9.51	8.6	89.1	
K5725-02	SO-2	51	1.16	9.66	8.44	85.6	
K5725-03	SO-3	52	1.11	9.88	8.33	82.3	
K5725-04	SO-4	53	1.11	9.72	8.38	84.4	
K5725-06	SO-1	66	1.15	9.51	8.6	89.1	
K5725-07	SO-2	67	1.16	9.66	8.44	85.6	
K5725-08	SO-3	68	1.11	9.88	8.33	82.3	
K5725-09	SO-4	69	1.11	9.72	8.38	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$