



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:10
In Date: 11/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 11/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106130

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5707-01	MH13Z4	1	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-02	MH13Z5	2	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-03	MH13Z6	3	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-04	MH13Z7	4	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-05	MH1400	5	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-06	MH1403	6	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-07	MH1404	7	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-08	MH1406	8	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-09	MH1407	9	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-10	MH1408	10	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-11	MH1409	11	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-12	MH1410	12	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-13	MH1411	13	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-14	MH1412	14	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-15	MH1412D	15	1.00	2.00	2.00	100.0	100% SOLIDS
K5707-16	MH1412S	16	1.00	2.00	2.00	100.0	100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$