



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/11/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:20
In Date: 11/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:22
Out Date: 11/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106188

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5772-01	C0AN4	1	1.14	9.93	8.86	87.8	
K5772-02	C0AN7	2	1.15	9.96	9.08	90.0	
K5772-03	C0AN8	3	1.15	9.56	8.53	87.8	
K5772-04	C0AN9	4	1.15	9.6	8.25	84.0	
K5772-05	C0AP0	5	1.14	9.54	8.31	85.4	
K5772-06	C0AP1	6	1.11	9.72	7.93	79.2	
K5772-07	C0AP2	7	1.18	9.85	8.54	84.9	
K5772-08	C0AP3	8	1.19	9.85	8.28	81.9	
K5772-09	C0AP4	9	1.11	9.98	8.63	84.8	
K5772-10	C0AP6	10	1.15	9.79	8.23	81.9	
K5772-11	C0AP7	11	1.19	9.61	8.67	88.8	
K5772-12	C0AP8	12	1.12	9.52	8.61	89.2	
K5772-13	C0AP9	13	1.18	9.71	7.93	79.1	
K5772-14	C0AQ0	14	1.18	9.8	8.97	90.4	
K5772-15	C0AQ1	15	1.18	9.63	8.63	88.2	
K5772-16	C0AQ2	16	1.19	9.89	8.44	83.3	
K5772-17	C0AQ3	17	1.15	9.64	8.8	90.1	
K5772-18	C0AQ4	18	1.17	9.91	8.43	83.1	
K5772-19	C0AQ4MS	19	1.17	9.91	8.43	83.1	
K5772-20	C0AQ4MSD	20	1.17	9.91	8.43	83.1	
K5772-21	C0AQ5	21	1.15	9.97	8.6	84.5	
K5772-22	C0AQ6	22	1.18	9.65	8.33	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$