



PERCENT SOLID

Supervisor: apatel
Analyst: Dodley
Date: 11/12/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:30
In Date: 11/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 11/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106195

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5776-01	C0AS8	1	1.15	9.85	8.75	87.4	
K5776-02	C0AS9	2	1.14	9.57	7.73	78.2	
K5776-03	C0AT0	3	1.14	9.75	8.8	89.0	
K5776-04	C0AT1	4	1.15	9.71	8.49	85.7	
K5776-05	C0AT2	5	1.15	9.83	7.3	70.9	
K5776-06	C0AT3	6	1.14	9.6	8.25	84.0	
K5776-07	C0AT4	7	1.15	9.81	8.24	81.9	
K5776-08	C0AT5	8	1.15	9.73	8.44	85.0	
K5776-09	C0AT6	9	1.15	9.62	8.05	81.5	
K5776-10	C0AT7	10	1.15	9.8	8.25	82.1	
K5776-11	C0AT8	11	1.15	9.61	8.76	90.0	
K5776-12	C0AT9	12	1.13	9.55	7.7	78.0	
K5776-13	C0AW0	13	1.13	9.67	8.04	80.9	
K5776-14	C0AW1	14	1.13	9.74	7.88	78.4	
K5776-15	C0AW2	15	1.14	9.74	8.47	85.2	
K5776-16	C0AW3	16	1.15	9.87	8.48	84.1	
K5776-17	C0AW4	17	1.13	9.83	8.36	83.1	
K5776-18	C0AW5	18	1.14	9.82	8.9	89.4	
K5776-19	C0AW7	19	1.14	9.78	8.43	84.4	
K5776-20	C0AX1	20	1.14	9.77	8.65	87.0	
K5776-21	C0AX1MS	21	1.14	9.77	8.65	87.0	
K5776-22	C0AX1MSD	22	1.14	9.77	8.65	87.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$