



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 11/14/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:26
In Date: 11/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 07:30
Out Date: 11/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106252

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5563-18	STOCKPILE	33	1.14	9.9	7.05	67.5	
K5563-19	ESW	34	1.14	9.8	8.17	81.2	
K5563-20	TP-1	35	1.14	9.56	7.37	74.0	
K5563-21	TP-2	36	1.13	9.85	7.54	73.5	
K5668-05	HD-02-111319	37	1.15	9.87	8.6	85.4	
K5668-06	HD-02-111319	38	1.14	9.7	8.48	85.7	
K5806-01	STOCK-PILE	1	1.13	9.75	8.57	86.3	
K5807-01	P001-CC01-191111-01	2	1.15	9.56	8.91	92.3	
K5807-03	P001-CC02-191111-01	3	1.16	9.96	9.16	90.9	
K5807-05	P001-CC03-191111-01	4	1.15	9.69	8.32	84.0	
K5807-07	P001-CC04-191111-01	5	1.15	9.65	9.06	93.1	
K5807-09	P001-CC05-191111-01	6	1.15	9.85	9.21	92.6	
K5807-11	P001-CC06-191111-01	7	1.16	9.81	8.82	88.6	
K5807-13	P001-CC07-191111-01	8	1.15	9.63	8.43	85.8	
K5807-15	P001-CC08-191112-01	9	1.14	9.58	8.94	92.4	
K5807-17	P001-CC09-191112-01	10	1.15	9.82	9.31	94.1	
K5807-19	P001-CC10-191112-01	11	1.16	9.74	8.95	90.8	
K5807-21	P001-CC11-191112-01	12	1.14	9.6	8.92	92.0	
K5807-23	P001-CC12-191112-01	13	1.15	9.94	9.32	92.9	
K5807-26	P001-CC13-191112-01	14	1.15	9.7	9.38	96.3	
K5808-01	P001-SB01-0006-191112-01	15	1.14	9.66	9.22	94.8	
K5808-02	P001-SB01-4248-191112-01	16	1.14	9.66	9.12	93.7	
K5808-03	P001-SB02-0006-191112-01	17	1.15	9.63	8.89	91.3	
K5808-04	P001-SB02-4248-191112-01	18	1.14	9.9	9.39	94.2	
K5808-05	P001-SB03-0006-191112-01	19	1.15	9.8	9.35	94.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 11/14/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:26
In Date: 11/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 07:30
Out Date: 11/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106252

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5808-06	P001-SB04-0006-191112-01	20	1.15	9.98	9.48	94.3	
K5808-07	P001-SB04-4248-191112-01	21	1.14	9.77	9.02	91.3	
K5808-08	P001-SB05-0006-191112-01	22	1.13	9.89	9.06	90.5	
K5808-09	P001-SB05-4248-191112-01	23	1.14	9.9	9.36	93.8	
K5808-10	P001-SB06-0006-191112-01	24	1.14	9.77	9.15	92.8	
K5808-11	P001-SB06-4248-191112-01	25	1.13	9.57	8.73	90.0	
K5808-12	K5808-11MS	26	1.13	9.57	8.73	90.0	
K5808-13	K5808-11MSD	27	1.13	9.57	8.73	90.0	
K5808-14	P001-SB06-4248-191112-02	28	1.13	9.89	9.33	93.6	
K5808-15	P001-SB07-0006-191112-01	29	1.13	9.86	9.51	96.0	
K5808-16	P001-SB07-4248-191112-01	30	1.15	9.77	9.39	95.6	
K5808-17	P001-SB08-0006-191111-01	31	1.14	9.8	9.37	95.0	
K5808-18	P001-SB08-4248-191111-01	32	1.16	9.88	9.33	93.7	
K5813-01	SP-B	39	1.15	9.76	7.33	71.8	
K5814-01	20688	40	1.00	2.00	2.00	100.0	oil sample 100% solids
K5818-01	30525	41	1.00	2.00	2.00	100.0	100% solids
K5818-03	40180	42	1.00	2.00	2.00	100.0	100% solids
K5826-01	TP-15-S1	43	1.15	9.98	8.94	88.2	
K5826-02	TP-15-S2	44	1.15	9.78	8.79	88.5	
K5826-03	TP-15-S	45	1.15	9.96	8.8	86.8	
K5826-05	TP-14B-S1	46	1.15	9.74	8.82	89.3	
K5826-06	TP-14B-S2	47	1.14	9.91	9.00	89.6	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 11/14/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:26
In Date: 11/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 07:30
Out Date: 11/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106252

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5826-07	TP-14B-S	48	1.14	9.58	8.67	89.2	
K5827-01	S-1	49	1.15	9.72	7.72	76.7	
K5827-02	S-2	50	1.14	9.7	7.35	72.5	
K5830-01	B-4- (2-4)	51	1.13	9.57	7.8	79.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$