



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 11/17/2019

OVENTEMP IN Celsius(°C): 104

Time IN: 12:00

In Date: 11/16/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102

Time OUT: 07:00

Out Date: 11/17/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB106332

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5879-01	MH1446	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-02	MH1447	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-03	MH1448	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-04	MH1449	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-05	MH1450	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-06	MH1451	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-07	MH1452	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-08	MH1453	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-09	MH1454	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-10	MH1455	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-11	MH1456	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-12	MH1457	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-13	MH1458	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-14	MH1459	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-15	MH1460	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-16	MH1461	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-17	MH1462	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-18	MH1463	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-19	MH1464	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-20	MH1466	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-21	MH1466D	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5879-22	MH1466S	22	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$