



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 11/17/2019

OVENTEMP IN Celsius(°C): 107

Time IN: 12:04

In Date: 11/16/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102

Time OUT: 07:05

Out Date: 11/17/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB106333

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5880-01	MH1465	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-02	MH1467	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-03	MH1468	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-04	MH1469	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-05	MH1470	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-06	MH1471	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-07	MH1472	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-08	MH1473	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-09	MH1474	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-10	MH1475	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-11	MH1476	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-12	MH1477	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-13	MH1478	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-14	MH1479	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-15	MH1480	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-16	MH1481	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-17	MH1482	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-18	MH1483	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-19	MH14A7	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-20	MH14A8	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-21	MH14A8D	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5880-22	MH14A8S	22	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$