



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 11/17/2019

OVENTEMP IN Celsius(°C): 107

Time IN: 12:09

In Date: 11/16/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 07:10

Out Date: 11/17/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB106334

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5881-01	MH1484	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-02	MH1485	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-03	MH1486	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-04	MH1487	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-05	MH1488	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-06	MH1489	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-07	MH1490	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-08	MH1491	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-09	MH1492	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-10	MH1493	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-11	MH1494	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-12	MH1495	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-13	MH1495D	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-14	MH1495S	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-15	MH1496	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-16	MH1497	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-17	MH1498	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-18	MH1499	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-19	MH14A0	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-20	MH14A1	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5881-21	MH14A2	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$