



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5669-01	CL-02-111519	1	1.15	9.96	9.29	92.4	
K5669-02	CL-02-111519	2	1.15	9.85	9.06	90.9	
K5669-03	CL-02-111519-1	3	1.15	9.64	9.28	95.8	
K5669-04	CL-02-111519-2	4	1.12	9.71	8.72	88.5	
K5669-05	CL-02-111519-3	5	1.15	9.95	9.3	92.6	
K5669-06	CL-02-111519-4	6	1.14	9.6	8.78	90.3	
K5669-07	CL-02-111519-5	7	1.14	9.94	9.16	91.1	
K5676-16	SU-01-111519	8	1.16	9.77	8.52	85.5	
K5676-17	SU-01-111519	9	1.12	9.93	9.33	93.2	
K5676-18	SU-01-111519-1	10	1.14	9.93	8.98	89.2	
K5676-19	SU-01-111519-2	11	1.14	9.93	9.26	92.4	
K5839-01	PHIPPS-BH-NW-01	12	1.15	9.94	9.24	92.0	
K5839-02	PHIPPS-BH-NW-01	13	1.16	9.73	8.94	90.8	
K5839-03	PHIPPS-BH-NW-02	14	1.19	9.62	8.91	91.6	
K5839-04	PHIPPS-BH-NW-02	15	1.11	9.54	8.99	93.5	
K5839-05	PHIPPS-BH-NW-03	16	1.18	9.7	9.45	97.1	
K5839-06	PHIPPS-BH-NW-03	17	1.16	9.87	9.61	97.0	
K5839-07	PHIPPS-BH-SW-01	18	1.17	9.89	9.26	92.8	
K5839-08	PHIPPS-BH-SW-01	19	1.18	9.97	9.33	92.7	
K5839-09	PHIPPS-BH-SW-02	20	1.16	9.74	8.66	87.4	
K5839-10	PHIPPS-BH-SW-02	21	1.13	9.94	9.64	96.6	
K5840-01	PHIPPS-BH-SW-03	22	1.18	9.91	9.6	96.4	
K5840-02	PHIPPS-BH-SW-03	23	1.15	9.5	9.06	94.7	
K5840-03	PHIPPS-BH-NE-01	24	1.19	9.88	9.06	90.6	
K5840-04	PHIPPS-BH-NE-01	25	1.16	9.67	9.22	94.7	
K5840-05	PHIPPS-BH-NE-02	26	1.18	9.75	8.43	84.6	
K5840-06	PHIPPS-BH-NE-02	27	1.15	9.97	9.04	89.5	
K5840-07	PHIPPS-BH-NE-03	28	1.15	9.98	8.85	87.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5840-08	PHIPPS-BH-NE-03	29	1.19	9.98	8.91	87.8	
K5840-09	PHIPPS-BH-SE-01	30	1.15	9.95	8.9	88.1	
K5840-10	PHIPPS-BH-SE-01	31	1.14	9.94	8.66	85.5	
K5840-11	PHIPPS-BH-SE-02	32	1.18	9.96	9.37	93.3	
K5840-12	PHIPPS-BH-SE-02	33	1.16	9.84	8.39	83.3	
K5840-13	PHIPPS-BH-SE-03	34	1.11	9.71	9.24	94.5	
K5840-14	PHIPPS-BH-SE-03	35	1.1	9.66	9.26	95.3	
K5841-01	P001-OSB01-0006-191113-01	36	1.12	9.75	8.92	90.4	
K5841-02	K5841-01MS	37	1.12	9.75	8.92	90.4	
K5841-04	P001-OSB01-0006-191113-02	38	1.14	9.7	9.51	97.8	
K5841-05	P001-OSB01-0612-191113-01	39	1.15	9.95	9.1	90.3	
K5841-06	P001-OSB01-4248-191113-01	40	1.13	9.69	9.00	91.9	
K5841-07	P001-OSB02-0006-191113-01	41	1.1	9.89	8.9	88.7	
K5841-08	P001-OSB02-0612-191113-01	42	1.15	9.96	9.36	93.2	
K5841-09	P001-OSB02-4248-191113-01	43	1.14	9.59	9.15	94.8	
K5841-10	P001-SB09-0006-191113-01	44	1.12	9.82	9.46	95.9	
K5841-11	P001-SB09-4248-191113-01	45	1.44	9.98	9.32	92.3	
K5841-12	P001-SB10-0006-191113-01	46	1.12	9.97	9.48	94.5	
K5841-13	P001-SB10-4248-191113-01	47	1.13	9.74	9.24	94.2	
K5841-14	P001-SB11-0006-191113-01	48	1.15	9.52	9.00	93.8	
K5841-15	P001-SB11-4248-191113-01	49	1.16	9.97	9.57	95.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5841-16	P001-SB12-0006-191113-01	50	1.19	9.69	9.33	95.8	
K5841-17	P001-SB13-0006-191113-01	51	1.11	9.7	9.46	97.2	
K5841-18	P001-SB13-4248-191113-01	52	1.13	9.82	9.51	96.4	
K5857-01	10327	53	1.00	2.00	2.00	100.0	100% SOLIDS
K5861-01	EP-1	54	1.13	9.5	8.46	87.6	
K5861-02	EP-2	55	1.16	9.86	8.93	89.3	
K5861-03	EP-3	56	1.15	9.96	9.4	93.6	
K5861-04	EP-4	57	1.16	9.79	9.22	93.4	
K5861-05	EP-5	58	1.16	9.51	8.99	93.8	
K5861-06	EP-6	59	1.13	9.5	9.16	95.9	
K5861-07	EP-7	60	1.19	9.65	9.12	93.7	
K5861-08	EP-8	61	1.19	9.82	9.55	96.9	
K5861-09	EP-9	62	1.19	9.97	9.48	94.4	
K5861-10	EP-10	63	1.11	9.83	9.47	95.9	
K5861-14	EP-11	64	1.15	9.98	9.3	92.3	
K5861-15	EP-12	65	1.14	9.94	8.98	89.1	
K5861-16	EP-13	66	1.15	9.95	9.12	90.6	
K5861-17	EP-14	67	1.14	9.74	9.18	93.5	
K5861-18	EP-15	68	1.14	9.97	9.76	97.6	
K5861-19	EP-16	69	1.16	9.84	9.44	95.4	
K5861-20	EP-17	70	1.15	9.96	9.46	94.3	
K5861-21	EP-18	71	1.11	9.52	9.1	95.0	
K5861-22	EP-19	72	1.15	9.97	9.38	93.3	
K5861-23	EP-20	73	1.12	9.9	9.53	95.8	
K5865-01	4-(2-4)	74	1.12	9.9	8.76	87.0	
K5865-02	4-(4-6)	75	1.14	9.78	9.14	92.6	
K5865-03	4-(6-8)	76	1.12	9.58	9.12	94.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5865-04	4-(8-10)	77	1.14	9.82	9.06	91.2	
K5865-05	4-(10-12)	78	1.13	9.81	9.32	94.4	
K5865-06	4-(12-14)	79	1.18	9.97	9.66	96.5	
K5865-07	4-(14-17)	80	1.13	9.77	9.36	95.3	
K5865-08	3-(2-4)	81	1.15	9.94	9.03	89.6	
K5865-09	3-(4-6)	82	1.18	9.91	9.02	89.8	
K5865-10	3-(6-8)	83	1.15	9.58	8.7	89.6	
K5865-11	3-(8-10)	84	1.14	9.8	8.98	90.5	
K5865-12	3-(10-12)	85	1.16	9.98	9.55	95.1	
K5865-13	3-(12-14)	86	1.15	9.96	9.67	96.7	
K5865-14	3-(14-16)	87	1.14	9.77	8.88	89.7	
K5865-15	3-(16-18)	88	1.13	9.98	9.44	93.9	
K5865-16	2-(2-4)	89	1.16	9.97	9.61	95.9	
K5865-17	2-(4-6)	90	1.13	9.68	8.92	91.1	
K5865-18	2-(6-8)	91	1.11	9.88	9.02	90.2	
K5865-19	2-(8-10)	92	1.13	9.72	9.31	95.2	
K5865-20	2-(10-12)	93	1.14	9.8	9.04	91.2	
K5865-21	2-(12-14)	94	1.15	9.68	9.06	92.7	
K5865-22	2-(14-16)	95	1.13	9.96	9.26	92.1	
K5865-23	2-(16-19)	96	1.14	9.78	9.13	92.5	
K5868-01	391107818	97	1.15	9.74	2.78	19.0	SLUDGE SAMPLE
K5869-01	391187713.1	98	1.16	9.86	2.48	15.2	SLUDGE SAMPLE
K5870-01	391107704	99	1.15	9.6	1.53	4.5	SLUDGE SAMPLE
K5871-01	391108035	100	1.17	9.98	1.49	3.6	SLUDGE SAMPLE
K5872-01	391112008	101	1.15	9.96	1.65	5.7	SLUDGE SAMPLE
K5876-01	1-7TH-FL-GYM-NE-CORNER	102	1.00	2.00	2.00	100.0	CLUCK SAMPLE 100% SOLIDS
K5876-02	2-7TH-FL-GYM-SE-CORNER	103	1.00	2.00	2.00	100.0	CLUCK SAMPLE 100% SOLIDS
K5893-01	COMP-RE	104	1.11	9.92	8.55	84.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5900-01	30695	105	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-03	30696	106	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-05	30697	107	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-07	30698	108	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-09	30699	109	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-11	30700	110	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-13	30724	111	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-15	30725	112	1.00	2.00	2.00	100.0	100% SOLIDS
K5900-17	30726	113	1.00	2.00	2.00	100.0	100% SOLIDS
K5903-01	EAST-END-EXTRA-SAMPLE	114	1.14	9.89	9.51	95.7	
K5904-01	1-(2-4)	115	1.14	9.71	9.21	94.2	
K5904-02	1-(4-6)	116	1.18	9.66	9.00	92.2	
K5904-03	1-(6-8)	117	1.13	9.62	9.03	93.1	
K5904-04	1-(8-10)	118	1.16	9.97	9.81	98.2	
K5904-05	1-(10-12)	119	1.1	9.88	9.5	95.7	
K5904-06	1-(12-14)	120	1.14	9.98	9.63	96.0	
K5904-07	1-(14-16)	121	1.14	9.79	9.41	95.6	
K5904-08	1-(16-18)	122	1.13	9.75	8.41	84.5	
K5904-09	1-(18-20)	123	1.15	9.92	9.04	90.0	
K5904-10	11-(0-2)	124	1.14	9.84	9.31	93.9	
K5904-11	11-(2-4)	125	1.44	9.63	8.2	82.5	
K5904-12	11-(4-7)	126	1.14	9.66	7.8	78.2	
K5904-13	12-(0-2)	127	1.16	9.61	8.68	89.0	
K5904-14	12-(2-4)	128	1.18	9.92	9.00	89.5	
K5904-15	12-(4-6)	129	1.18	9.75	7.75	76.7	
K5904-16	13-(0-2)	130	1.2	9.98	9.12	90.2	
K5904-17	13-(2-5)	131	1.16	9.64	8.1	81.8	
K5904-18	14-(0-2)	132	1.16	9.72	8.87	90.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5904-19	14-(2-4)	133	1.18	9.97	9.06	89.6	
K5905-01	RT-3037	134	1.15	9.95	8.28	81.0	
K5905-02	RT-3037	135	1.15	9.95	8.28	81.0	
K5906-01	NA-SP	136	1.18	9.96	9.33	92.8	
K5906-02	NA-SP	137	1.18	9.96	9.33	92.8	
K5906-03	NA-SP	138	1.18	9.96	9.33	92.8	
K5907-01	99K-1B	139	1.12	9.93	8.32	81.7	
K5907-02	99K-1C	140	1.15	9.94	8.13	79.4	
K5907-03	99K-2B	141	1.14	9.81	8.2	81.4	
K5907-04	99K-2C	142	1.1	9.51	7.62	77.5	
K5907-05	99K-3A	143	1.13	9.78	8.09	80.5	
K5907-06	99K-3B	144	1.19	9.83	7.32	70.9	
K5907-07	99K-3C	145	1.16	9.82	7.95	78.4	
K5907-08	99K-4C	146	1.19	9.79	7.81	77.0	
K5907-09	99K-5C	147	1.14	9.5	7.6	77.3	
K5907-10	99K-6A	148	1.19	9.58	8.03	81.5	
K5907-11	99K-6B	149	1.11	9.55	7.64	77.4	
K5907-12	99K-6C	150	1.16	9.66	8.08	81.4	
K5907-13	99K-7B	151	1.18	9.66	7.93	79.6	
K5907-14	99K-7C	152	1.16	9.77	7.98	79.2	
K5907-15	99K-8B	153	1.18	9.98	8.6	84.3	
K5907-16	99K-8C	154	1.18	9.97	8.47	82.9	
K5907-17	99K-DUP-1	155	1.18	9.76	8.06	80.2	
K5907-18	99K-DUP-2	156	1.15	9.95	8.43	82.7	
K5907-19	99K-DUP-3	157	1.12	9.62	7.81	78.7	
K5908-01	B-8D-(12-18)	158	1.18	9.77	8.56	85.9	
K5908-02	B-8-1-(6-12)	159	1.18	9.79	8.94	90.1	
K5908-03	B-8-2-(6-12)	160	1.18	9.51	8.73	90.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/18/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 11/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:15
Out Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5908-04	B-8-3-(6-12)	161	1.12	9.92	9.24	92.3	
K5909-05	SVOC-GPC-BLANK	162	1.00	2.00	2.00	100.0	
K5909-06	PEST-GPC-BLANK	163	1.00	2.00	2.00	100.0	
K5909-07	PEST-GPC-BLANK-SPIKE	164	1.00	2.00	2.00	100.0	
K5909-08	PCB-GPC-BLANK	165	1.00	2.00	2.00	100.0	
K5909-09	PCB-GPC-BLANK-SPIKE	166	1.00	2.00	2.00	100.0	
K5909-10	SVOC-GPC2-BLANK	167	1.00	2.00	2.00	100.0	
K5909-11	PEST-GPC2-BLANK	168	1.00	2.00	2.00	100.0	
K5909-12	PEST-GPC2-BLANK-SPIKE	169	1.00	2.00	2.00	100.0	
K5909-13	PCB-GPC2-BLANK	170	1.00	2.00	2.00	100.0	
K5909-14	PCB-GPC2-BLANK -SPIKE	171	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$