



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 11/19/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 11/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106357

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5894-01	MBFQJ5	1	1.14	9.91	9.77	98.4	
K5894-02	MBFQJ5D	2	1.14	9.91	9.77	98.4	
K5894-03	MBFQJ5S	3	1.14	9.91	9.77	98.4	
K5894-04	MBFQJ6	4	1.15	9.82	9.53	96.7	
K5894-05	MBFQJ7	5	1.13	9.56	9.3	96.9	
K5894-06	MBFQJ8	6	1.15	9.94	9.78	98.2	
K5894-07	MBFQJ9	7	1.18	9.67	9.48	97.8	
K5894-08	MBFQK0	8	1.19	9.61	9.38	97.3	
K5894-09	MBFQK1	9	1.18	9.61	7.38	73.5	
K5894-10	MBFQK2	10	1.17	9.58	7.37	73.7	
K5894-11	MBFQK3	11	1.15	9.7	7.74	77.1	
K5894-12	MBFQK4	12	1.14	9.98	9.02	89.1	
K5894-13	MBFQK5	13	1.17	9.66	8.77	89.5	
K5894-14	MBFQK6	14	1.11	9.6	8.68	89.2	
K5894-15	MBFQK7	15	1.15	9.58	7.79	78.8	
K5894-16	MBFQK8	16	1.14	9.67	7.84	78.5	
K5894-17	MBFQK9	17	1.15	9.68	7.92	79.4	
K5894-18	MBFQL0	18	1.17	9.69	8.16	82.0	
K5894-19	MBFQL2	19	1.15	9.95	8.33	81.6	
K5894-20	MBFQL3	20	1.11	9.66	9.19	94.5	
K5894-21	MBFQL4	21	1.14	9.62	9.14	94.3	
K5894-22	MBFQL5	22	1.13	9.52	9.12	95.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$