



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 11/19/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 11/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:47
Out Date: 11/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106362

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5913-01	MJLNT4	1	1.15	9.94	7.44	71.6	
K5913-02	MJLNT5	2	1.14	9.98	7.51	72.1	
K5913-03	MJLNT6	3	1.14	9.96	7.06	67.1	
K5913-04	MJLNT7	4	1.11	9.88	6.86	65.6	
K5913-05	MJLNT7D	5	1.11	9.88	6.86	65.6	
K5913-06	MJLNT7S	6	1.11	9.88	6.86	65.6	
K5913-07	MJLNT8	7	1.15	9.59	8.16	83.1	
K5913-08	MJLNT9	8	1.15	9.96	6.26	58.0	
K5913-09	MJLNW0	9	1.19	9.71	5.97	56.1	
K5913-10	MJLNW1	10	1.14	9.83	6.77	64.8	
K5913-11	MJLNW2	11	1.15	9.96	7.51	72.2	
K5913-12	MJLNW3	12	1.18	9.96	7.68	74.0	
K5913-13	MJLNW4	13	1.13	9.9	7.26	69.9	
K5913-14	MJLNW5	14	1.1	9.76	7.42	73.0	
K5913-15	MJLNW6	15	1.18	9.66	6.16	58.7	
K5913-16	MJLNW7	16	1.13	9.57	7.15	71.3	
K5913-17	MJLNW8	17	1.18	9.88	6.78	64.4	
K5913-18	MJLNW9	18	1.15	9.93	6.48	60.7	
K5913-19	MJLNX0	19	1.13	9.97	5.88	53.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$