

LB106415

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AK

Instrument ID : Konelab

11/20/2019 16:33

Test: ~~Total~~ CN

Reactive
ug/L

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	97.168	0.0	0.102	
ICB1	0.679	0.0	0.002	
CCV1	251.291	0.0	0.261	
CCB1	0.541	0.0	0.002	
PB124953BL	0.397	0.0	0.002	
K5916-02	0.450	0.0	0.002	
K5916-02DUP	0.423	0.0	0.002	
K5933-05	0.537	0.0	0.002	
K5933-11	0.444	0.0	0.002	
K5934-01	0.121	0.0	0.002	
K5934-02	0.434	0.0	0.002	
K5938-02	0.093	0.0	0.002	
CCV2	266.481	0.0	0.277	
CCB2	0.395	0.0	0.002	

N	14
Mean	44.247
SD	94.5411
CV%	213.67

Aquakem v. 7.2AQ1

AK

Results from time period:

Wed Nov 20 16:23:32 2019

Wed Nov 20 16:27:25 2019

Sample Id	Sam/Ctr/c/	Test short Test type	Result	Result unit	Result date and time	Stat
0.0PPBCN	A	Reactive Cl P	0.0223	µg/l	11/20/2019 14:57:18	
5.0PPBCN	A	Reactive Cl P	5.6674	µg/l	11/20/2019 14:57:19	
10PPBCN	A	Reactive Cl P	11.4803	µg/l	11/20/2019 14:57:20	
50PPBCN	A	Reactive Cl P	46.7606	µg/l	11/20/2019 14:57:21	
100PPBCN	A	Reactive Cl P	96.7582	µg/l	11/20/2019 14:57:22	
250PPBCN	A	Reactive Cl P	256.7502	µg/l	11/20/2019 14:57:23	
500PPBCN	A	Reactive Cl P	497.5609	µg/l	11/20/2019 14:57:24	
ICV1	S	Reactive Cl P	97.1676	µg/l	11/20/2019 16:23:32	
ICB1	S	Reactive Cl P	0.6789	µg/l	11/20/2019 16:23:33	
CCV1	S	Reactive Cl P	251.2914	µg/l	11/20/2019 16:23:34	
CCB1	S	Reactive Cl P	0.5412	µg/l	11/20/2019 16:23:35	
PB124953BL	S	Reactive Cl P	0.3971	µg/l	11/20/2019 16:23:36	
K5916-02	S	Reactive Cl P	0.4497	µg/l	11/20/2019 16:23:37	
K5916-02DUP	S	Reactive Cl P	0.4232	µg/l	11/20/2019 16:23:38	
K5933-05	S	Reactive Cl P	0.5373	µg/l	11/20/2019 16:23:39	
K5933-11	S	Reactive Cl P	0.4437	µg/l	11/20/2019 16:23:40	
K5934-01	S	Reactive Cl P	0.1206	µg/l	11/20/2019 16:23:41	
K5934-02	S	Reactive Cl P	0.4336	µg/l	11/20/2019 16:23:42	
K5938-02	S	Reactive Cl P	0.0932	µg/l	11/20/2019 16:27:23	
CCV2	S	Reactive Cl P	266.4807	µg/l	11/20/2019 16:27:24	
CCB2	S	Reactive Cl P	0.3951	µg/l	11/20/2019 16:27:25	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AK

Instrument ID : Konelab

11/20/2019 15:04


Test Total CN

Accepted 11/20/2019 15:04

Factor 966.5
Bias 0.001

Coeff. of det. 0.999634

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.001	0.0223	0.0000	
2	5.0PPBCN	0.007	5.6674	5.0000	
3	10PPBCN	0.013	11.4803	10.0000	
4	50PPBCN	0.050	46.7606	50.0000	
5	100PPBCN	0.102	96.7582	100.0000	
6	250PPBCN	0.267	256.7502	250.0000	
7	500PPBCN	0.516	497.5609	500.0000	