



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 11/25/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 11/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106451

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5990-01	MH14H9	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-02	MH14J0	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-03	MH14J1	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-04	MH14J2	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-05	MH14J3	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-06	MH14J4	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-07	MH14J4D	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-08	MH14J4S	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-09	MH14J5	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-10	MH14J6	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-11	MH14J7	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-12	MH14J8	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-13	MH14J9	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-14	MH14K0	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-15	MH14K1	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-16	MH14K2	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-17	MH14K3	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-18	MH14K8	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-19	MH14K9	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-20	MH14L0	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-21	MH14L1	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5990-22	MH14L2	22	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$