



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/5/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 12/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 09:13
Out Date: 12/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106587

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5956-04	FAIRMOUNT-FOUNDATION-3	1	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K6102-03	391126164	2	1.14	9.97	1.16	0.2	SLUDGE SAMPLE
K6103-01	391126119	3	1.15	9.96	1.2	0.6	SLUDGE SAMPLE
K6112-01	OK-02-120419	16	1.13	9.61	9.39	97.4	
K6112-02	OK-02-120419	17	1.15	9.9	9.43	94.6	
K6112-03	OK-02-120419-1	18	1.14	9.63	9.1	93.8	
K6119-01	SB-3- (0-2)	4	1.14	9.9	9.39	94.2	
K6119-02	SB-3- (8-10)	5	1.15	9.92	9.5	95.2	
K6119-03	SB-6- (1-3)	6	1.15	9.9	9.23	92.3	
K6119-04	SB-6- (7-9)	7	1.14	9.87	9.35	94.0	
K6119-05	SB-4- (0-2)	8	1.13	9.69	8.74	88.9	
K6119-06	SB-4- (8-10)	9	1.14	9.56	9.11	94.7	
K6119-07	SB-5- (0-2)	10	1.14	9.94	8.63	85.1	
K6119-08	SB-5- (8-10)	11	1.19	9.68	9.07	92.8	
K6119-09	SB-2- (0-2)	12	1.18	9.92	8.79	87.1	
K6119-10	SB-2- (8-10)	13	1.12	9.61	8.74	89.8	
K6119-11	SB-1- (0-2)	14	1.17	9.6	9.03	93.2	
K6119-12	SB-1- (8-10)	15	1.18	9.98	9.33	92.6	
K6133-01	ROLL-OFF-362	19	1.16	9.51	8.52	88.1	
K6135-01	VOC-34	20	1.16	6.00	5.56	90.9	
K6135-02	COMP-34	21	1.15	9.81	9.26	93.6	
K6135-03	VOC-35	22	1.15	6.48	6.03	91.6	
K6135-04	COMP-35	23	1.15	9.53	8.81	91.4	
K6135-05	VOC-36	24	1.16	6.25	5.69	89.0	
K6135-06	COMP-36	25	1.16	9.57	8.81	91.0	
K6141-01	CMJ-01-1.5-2.0	26	1.14	9.61	8.53	87.2	
K6141-02	CMJ-01-2.5-3.0	27	1.13	9.54	8.4	86.4	
K6141-03	CMJ-01-7.5-8.0	28	1.15	9.89	8.6	85.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/5/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 12/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 09:13
Out Date: 12/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106587

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6141-04	CMJ-01-COMP	29	1.15	9.78	8.78	88.4	
K6142-01	1D	30	1.13	9.54	7.88	80.3	
K6142-04	2D	31	1.14	9.94	8.43	82.8	
K6142-07	3D	32	1.15	9.96	8.45	82.9	
K6142-10	4D	33	1.14	9.76	8.28	82.8	
K6142-11	4D-DUP	34	1.14	9.77	8.01	79.6	
K6142-14	5D	35	1.14	9.85	7.96	78.3	
K6142-17	6D	36	1.15	9.81	8.14	80.7	
K6142-20	7D	37	1.12	9.96	8.21	80.2	
K6143-04	8D	38	1.15	9.58	6.87	67.9	
K6143-07	9D	39	1.16	9.83	7.45	72.5	
K6143-10	10D	40	1.15	9.79	8.05	79.9	
K6143-14	11D	41	1.14	9.9	8.06	79.0	
K6143-17	12D	42	1.16	9.82	8.28	82.2	
K6143-20	13D	43	1.16	9.63	8.13	82.3	
K6144-03	14D	44	1.15	9.76	7.24	70.7	
K6144-04	14D-DUP	45	1.15	9.54	7.59	76.8	
K6144-07	15D	46	1.14	9.84	8.13	80.3	
K6144-10	16D	47	1.13	9.96	8.00	77.8	
K6145-01	TH-2	48	1.15	9.6	8.3	84.6	
K6145-03	TH-3	49	1.16	9.61	8.3	84.5	
K6145-05	TH-8	50	1.16	9.89	7.95	77.8	
K6147-01	P001-WC001	51	1.14	9.86	8.76	87.4	
K6147-02	P001-WC002	52	1.16	9.8	9.15	92.5	
K6147-03	P001-WC003	53	1.14	9.88	9.24	92.7	
K6147-04	P001-WC004	54	1.12	9.98	9.25	91.8	
K6147-05	P001-WC005	55	1.15	9.87	9.1	91.2	
K6147-06	P001-WC006	56	1.12	9.88	9.16	91.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/5/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 12/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 09:13
Out Date: 12/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106587

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6147-07	P001-WC007	57	1.15	9.51	8.67	90.0	
K6147-08	P001-WC008	58	1.15	9.87	9.03	90.4	
K6147-09	P001-WC009	59	1.13	9.6	8.74	89.8	
K6147-10	P001-WC010	60	1.15	9.55	8.64	89.2	
K6147-11	P001-WC011	61	1.15	9.62	8.56	87.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$