



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 12/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 12/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106601

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K6153-01	P001-CA42-0612-01	1	1.14	9.94	9.24	92.0	
K6153-02	P001-CA43-0612-01	2	1.18	9.96	9.27	92.1	
K6153-03	P001-CA44-0006-01	3	1.15	9.77	9.02	91.3	
K6153-04	P001-CB41-0612-01	4	1.15	9.71	9.09	92.8	
K6153-05	P001-CB42-0612-01	5	1.12	9.69	9.06	92.6	
K6153-06	P001-CB43-0612-01	6	1.11	9.52	9.17	95.8	
K6153-07	P001-CB44-0612-01	7	1.18	9.68	9.41	96.8	
K6153-08	P001-CB45-0612-01	8	1.16	9.83	8.98	90.2	
K6153-09	P001-CC40-0006-01	9	1.17	9.82	8.8	88.2	
K6153-10	P001-CC41-0612-01	10	1.13	9.55	8.7	89.9	
K6153-11	P001-CC42-1218-01	11	1.14	9.62	8.98	92.5	
K6153-12	P001-CC43-0612-01	12	1.14	9.75	9.44	96.4	
K6153-13	P001-CC44-1218-01	13	1.14	9.94	9.61	96.3	
K6153-14	P001-CC45-1218-01	14	1.18	9.62	9.05	93.2	
K6153-15	P001-CC46-0612-01	15	1.18	9.97	9.62	96.0	
K6153-16	P001-CC47-1218-01	16	1.14	9.83	9.16	92.3	
K6153-17	P001-CD40-0006-01	17	1.15	9.9	9.12	91.1	
K6153-18	P001-CD41-0006-01	18	1.19	9.7	8.82	89.7	
K6153-19	P001-CD42-1218-01	19	1.14	9.77	9.15	92.8	
K6153-20	P001-CE48-1218-01	20	1.15	9.77	9.14	92.7	
K6153-21	P001-CE48-1218-01D	21	1.15	9.77	9.14	92.7	
K6153-22	P001-CE48-1218-01S	22	1.15	9.77	9.14	92.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$