



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 12/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 12/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106914

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K6413-01	K005-AA-SB1	1	1.13	7.06	6.48	90.2	
K6413-02	K005-AA-SB2	2	1.15	7.53	6.88	89.8	
K6413-03	K005-AA-SB3	3	1.13	6.83	6.42	92.8	
K6413-04	K005-AA-SB4	4	1.14	6.06	5.22	82.9	
K6413-05	K005-AA-SB5	5	1.12	7.24	6.66	90.5	
K6413-06	K005-AA-WC-1	6	1.18	9.79	8.79	88.4	
K6413-07	K005-AA-WC-2	7	1.11	9.65	8.9	91.2	
K6413-08	K005-AB-SB1	8	1.15	5.95	5.48	90.2	
K6413-09	K005-AB-SB2	9	1.16	5.52	5.15	91.5	
K6413-10	K005-AB-SB3	10	1.13	6.05	5.69	92.7	
K6413-11	K005-AB-SB4	11	1.16	6.99	6.48	91.3	
K6413-12	K005-AB-SB5	12	1.12	6.19	5.82	92.7	
K6413-13	K005-AB-WC-1	13	1.13	9.81	8.94	90.0	
K6414-01	K183-SB-01	14	1.17	6.73	6.09	88.5	
K6414-02	K183-SB-02	15	1.18	7.2	6.56	89.4	
K6414-03	K183-SB-03	16	1.11	7.42	6.88	91.4	
K6414-04	K183-SB-04	17	1.18	7.62	6.67	85.2	
K6414-05	K183-SB-05	18	1.15	7.08	6.3	86.8	
K6414-06	K183-SB-06	19	1.17	5.55	5.42	97.0	
K6414-07	K183-WC-01	20	1.16	9.97	9.00	89.0	
K6414-08	K183-WC-02	21	1.14	9.6	8.81	90.7	
K6426-01	391218026.0	22	1.15	9.96	1.51	4.1	SLUDGE SAMPLE
K6427-01	391219032.0	23	1.18	9.82	2.81	18.9	SLUDGE SAMPLE
K6428-01	391218116.0	24	1.17	9.98	1.44	3.1	SLUDGE SAMPLE
K6429-01	391217172.0	25	1.19	9.97	1.3	1.3	SLUDGE SAMPLE
K6430-01	391218726.0	26	1.14	9.86	3.09	22.4	SLUDGE SAMPLE
K6431-01	SB-5- (0-2)	57	1.14	9.95	9.01	89.3	
K6431-02	SB-5- (4-5)	58	1.14	9.82	8.84	88.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 12/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 12/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106914

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6431-03	SB-4- (0-2)	59	1.13	9.61	8.44	86.2	
K6431-04	SB-4- (5-6)	60	1.16	9.7	8.33	84.0	
K6431-06	SB-3- (0-2)	61	1.13	9.89	8.36	82.5	
K6431-07	SB-3- (4-5)	62	1.1	9.78	8.83	89.1	
K6431-08	SB-2- (0-2)	63	1.11	9.8	8.33	83.1	
K6431-09	SB-2- (4-5)	64	1.15	9.84	8.92	89.4	
K6431-10	K6431-09MS	65	1.15	9.84	8.92	89.4	
K6431-11	K6431-09MSD	66	1.15	9.84	8.92	89.4	
K6431-12	SB-1- (0-2)	67	1.13	9.78	8.65	86.9	
K6431-13	SB-1- (4-5)	68	1.18	9.97	9.1	90.1	
K6431-14	SOIL-DUPLICATE	69	1.12	9.55	8.18	83.7	
K6433-01	01-A-01-B-01-C	27	1.00	2.00	2.00	100.0	CAULKING MATERIAL 100% SOLIDS
K6435-01	SB-4D- (21.5-23.5)	70	1.15	9.97	9.22	91.5	
K6435-02	SB-4E- (21.5-23.5)	71	1.17	9.94	8.43	82.8	
K6435-03	SB-3E- (09.5)	72	1.17	9.59	7.98	80.9	
K6435-04	SB-3E- (09.5-11.5)	73	1.15	9.88	8.13	80.0	
K6436-01	K009-SB-1	74	1.14	4.45	4.29	95.2	
K6436-02	K009-SB-2	75	1.13	5.92	5.59	93.1	
K6436-03	K009-SB-3	76	1.15	6.16	5.91	95.0	
K6436-04	K009-SB-4	77	1.15	6.64	6.21	92.2	
K6436-05	K009-SB-5	78	1.11	6.77	6.32	92.0	
K6436-06	K009-WC-1	79	1.15	9.56	8.8	91.0	
K6436-07	K009-WC-2	80	1.14	9.9	9.14	91.3	
K6437-01	GP-2- (1-3)	81	1.18	9.91	9.18	91.6	
K6437-02	GP-2- (3-5)	82	1.14	9.96	9.00	89.1	
K6437-03	GP-1- (1-3)	83	1.15	9.69	8.36	84.4	
K6437-04	GP-1- (3-5)	84	1.1	9.93	8.98	89.2	
K6437-05	GP-3- (1-3)	85	1.18	9.7	8.82	89.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 12/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 12/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106914

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6437-06	GP-3-(3-5)	86	1.14	9.83	8.59	85.7	
K6437-07	GP-5-(5-6)	87	1.15	9.96	8.68	85.5	
K6437-08	GP-5-(6-8)	88	1.18	9.71	7.44	73.4	
K6438-01	DUNR-BF001-01	28	1.14	9.58	9.1	94.3	
K6438-02	K6438-01MS	29	1.14	9.58	9.1	94.3	
K6438-03	K6438-01MSD	30	1.14	9.58	9.1	94.3	
K6438-04	DUNR-BF001-02	31	1.14	9.98	9.55	95.1	
K6438-05	DUNR-BF002-01	32	1.14	9.6	9.13	94.4	
K6438-06	DUNR-BF003-01	33	1.14	9.96	9.4	93.7	
K6438-07	DUNR-BF004-01	34	1.15	9.9	9.56	96.1	
K6438-08	DUNR-BF005-01	35	1.13	9.95	9.48	94.7	
K6438-09	DUNR-BF006-01	36	1.14	9.66	9.22	94.8	
K6438-10	DUNR-BF007-01	37	1.15	9.97	9.48	94.4	
K6438-11	DUNR-BF008-01	38	1.14	9.59	9.32	96.8	
K6438-12	DUNR-BF009-01	39	1.14	9.62	9.06	93.4	
K6439-01	LIED-BF001-01	40	1.14	9.67	8.99	92.0	
K6439-02	K6439-01MS	41	1.14	9.67	8.99	92.0	
K6439-03	K6439-01MSD	42	1.14	9.67	8.99	92.0	
K6439-04	LIED-BF001-02	43	1.13	9.95	9.43	94.1	
K6439-05	LIED-BF002-01	44	1.15	9.76	9.03	91.5	
K6439-06	LIED-BF003-01	45	1.15	9.64	8.85	90.7	
K6439-07	LIED-BF004-01	46	1.15	9.7	9.05	92.4	
K6439-08	LIED-BF005-01	47	1.12	9.65	9.11	93.7	
K6439-09	LIED-BF006-01	48	1.13	9.82	9.21	93.0	
K6439-10	LIED-BF007-01	49	1.14	9.98	9.06	89.6	
K6439-11	LIED-BF008-01	50	1.14	9.53	8.85	91.9	
K6439-12	LIED-BF009-01	51	1.14	9.63	8.78	90.0	
K6439-13	LIED-BF010-01	52	1.13	9.97	9.63	96.2	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 12/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 12/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106914

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6439-14	LIED-BF011-01	53	1.13	9.65	9.15	94.1	
K6439-15	LIED-BF012-01	54	1.14	9.92	9.61	96.5	
K6439-16	LIED-BF013-01	55	1.14	9.97	9.52	94.9	
K6439-17	LIED-BF014-01	56	1.15	9.76	9.48	96.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$