



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/31/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 12/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 12/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106963

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6462-01	K500-1	1	1.14	9.95	8.87	87.7	
K6462-02	K500-2	2	1.12	9.53	8.72	90.4	
K6462-03	K500-3	3	1.15	9.66	8.83	90.2	
K6462-04	K500-4	4	1.17	9.62	8.9	91.5	
K6462-05	K500-5	5	1.18	9.64	8.88	91.0	
K6463-01	GP-6- (1-3)	6	1.15	9.6	8.64	88.6	
K6463-02	GP-6- (4-6)	7	1.13	9.94	8.61	84.9	
K6463-03	GP-7- (1-3)	8	1.13	9.56	8.77	90.6	
K6463-04	GP-7- (4-6)	9	1.17	9.9	9.37	93.9	
K6463-05	GP-5C- (5-6)	10	1.16	9.69	8.4	84.9	
K6464-01	SB-4- (5-6)	11	1.15	9.96	8.13	79.2	
K6464-02	SB-5- (8-10)	12	1.18	9.85	8.6	85.6	
K6465-01	GP-11- (1-3)	13	1.18	9.58	8.52	87.4	
K6465-02	GP-11- (5-6)	14	1.13	9.83	8.37	83.2	
K6465-03	GP-10- (1-3)	15	1.17	9.67	8.61	87.5	
K6465-04	GP-10- (10-12)	16	1.17	9.68	7.07	69.3	
K6465-05	GP-9- (1-3)	17	1.1	9.88	8.87	88.5	
K6465-06	GP-9- (8-10)	18	1.19	9.85	8.19	80.8	
K6465-07	GP-4- (1-3)	19	1.11	9.77	8.74	88.1	
K6465-08	GP-4- (3-5)	20	1.18	9.94	8.68	85.6	
K6465-09	SS-8	21	1.12	9.74	8.47	85.3	
K6469-01	174-NC-TOPSOIL-001	22	1.13	9.76	6.37	60.7	
K6471-01	MC-123019-0-1-GRAB-B5	23	1.18	9.65	8.25	83.5	
K6471-02	MC-123019-0-2-COMP-B5	24	1.13	9.69	8.22	82.8	
K6471-03	MC-123019-4-5-GRAB-B5	25	1.14	9.51	8.16	83.9	
K6471-04	MC-123019-2-5-COMP-B5	26	1.13	9.94	8.62	85.0	
K6471-05	MC-123019-0-1-GRAB-B8	27	1.16	9.55	7.96	81.0	
K6471-06	MC-123019-0-2-COMP-B8	28	1.11	9.68	8.11	81.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/31/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 12/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 12/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106963

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6471-07	MC-123019-4-5-GRAB-B8	29	1.18	9.72	8.23	82.6	
K6471-08	MC-123019-2-5-COMP-B8	30	1.19	9.7	8.07	80.8	
K6471-09	MC-123019-0-1-GRAB-B2	31	1.18	9.75	8.38	84.0	
K6471-10	MC-123019-0-2-COMP-B2	32	1.11	9.94	8.7	86.0	
K6471-11	MC-123019-4-5-GRAB-B2	33	1.18	9.92	8.62	85.1	
K6471-12	MC-123019-2-5-COMP-B2	34	1.14	9.98	9.33	92.6	
K6471-13	MC-123019-0-1-GRAB-B1	35	1.14	9.68	8.79	89.6	
K6471-14	MC-123019-0-2-COMP-B1	36	1.17	9.69	8.8	89.6	
K6471-15	MC-123019-2-5-COMP-B1	37	1.1	9.84	8.75	87.5	
K6471-16	MC-123019-4-5-GRAB-B1	38	1.13	9.57	8.43	86.5	
K6472-01	K309-SB-01	39	1.13	7.54	7.04	92.2	
K6472-02	K309-SB-02	40	1.15	6.88	6.53	93.9	
K6472-03	K309-SB-03	41	1.18	7.17	6.3	85.5	
K6472-05	K309-SB-05	42	1.14	7.41	6.45	84.7	
K6472-06	K309-WC-01	43	1.14	9.58	8.6	88.4	
K6474-01	SB-01	44	1.18	9.78	8.97	90.6	
K6474-03	SB-02	45	1.18	9.94	9.11	90.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$