



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/31/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 12/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 12/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB106982

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6445-13	X2K09	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-14	X2K09D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-15	X2K09S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-16	X2K10	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-17	X2K10D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-18	X2K10S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-19	X2K11	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-20	X2K11D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-21	X2K11S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-22	X2K12	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-23	X2K12D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
K6445-24	X2K12S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$