



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 1/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 01/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107005

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6484-16	40096	1	1.12	9.92	7.22	69.3	
K6484-17	40097	2	1.18	9.67	6.86	66.9	
K6484-18	40099	3	1.18	9.62	5.82	55.0	
K6484-19	40100	4	1.15	9.95	7.51	72.3	
K6484-20	40101	5	1.12	9.71	6.1	58.0	
K6484-21	40102	6	1.17	9.58	5.6	52.7	
K6484-22	40103	7	1.13	7.3	3.52	38.7	
K6484-23	40104	8	1.13	9.88	6.29	59.0	
K6484-24	40105	9	1.16	9.72	4.26	36.2	
K6484-25	40106	10	1.15	9.9	5.93	54.6	
K6484-26	40107	11	1.11	9.92	5.57	50.6	
K6484-27	40108	12	1.17	9.98	7.19	68.3	
K6484-28	40109	13	1.17	9.97	7.08	67.2	
K6484-29	40110	14	1.17	9.98	7.53	72.2	
K6484-30	40111	15	1.15	9.79	6.61	63.2	
K6484-31	40112	16	1.16	9.96	7.46	71.6	
K6484-32	40113	17	1.16	9.92	6.25	58.1	
K6484-33	40114	18	1.11	9.73	5.00	45.1	
K6484-34	40115	19	1.16	9.59	5.93	56.6	
K6484-35	40116	20	1.15	6.3	3.79	51.3	
K6484-36	40117	21	1.14	9.91	7.43	71.7	
K6484-37	40118	22	1.19	9.98	6.08	55.6	
K6484-38	40119	23	1.13	9.7	4.49	39.2	
K6484-39	40120	24	1.18	9.64	6.9	67.6	
K6484-44	40096-40103	25	1.13	9.98	6.85	64.6	
K6484-46	40104-40111	26	1.13	9.8	6.12	57.6	
K6484-48	40112-40120	27	1.19	9.97	6.36	58.9	
K6484-50	30544-46-66	28	1.00	2.00	2.00	100.0	100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 1/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 01/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107005

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K6484-52	30544	29	1.00	2.00	2.00	100.0	100% SOLIDS
K6484-53	30546	30	1.00	2.00	2.00	100.0	100% SOLIDS
K6484-54	30566	31	1.11	9.97	7.74	74.8	
L1001-01	PAMT-R0R0-CORE-3-12-30-19	32	1.13	9.59	5.23	48.5	
L1001-02	RAMT-R0R0-COMP-1-12-30-19	33	1.13	9.85	4.73	41.3	
L1001-03	SC-CORE-3-12-31-19	34	1.16	9.96	4.11	33.5	
L1001-04	SC-1-12-31-19	35	1.13	9.73	5.15	46.7	
L1001-05	PAMT-COMP-1-12-31-19	36	1.14	9.96	6.09	56.1	
L1001-06	L1001-05MS	37	1.14	9.96	6.09	56.1	
L1001-07	L1001-05MSD	38	1.14	9.96	6.09	56.1	
L1001-08	PAMT-COMP-2-12-31-19	39	1.11	9.55	5.21	48.6	
L1008-01	CL-02-010220	48	1.14	9.89	9.18	91.9	
L1008-03	CL-02-010220	49	1.15	9.93	8.71	86.1	
L1008-04	CL-02-010220-1	50	1.14	9.88	9.00	89.9	
L1008-05	CL-02-010220-2	51	1.18	9.95	8.97	88.8	
L1008-06	CL-02-010220-3	52	1.12	9.58	9.46	98.6	
L1008-07	CL-02-010220-4	53	1.17	9.6	9.07	93.7	
L1009-01	NB-08-010220	54	1.17	9.91	8.84	87.8	
L1009-03	NB-08-010220	55	1.18	9.92	8.43	83.0	
L1009-04	NB-08-010220-1	56	1.13	9.94	8.3	81.4	
L1009-05	NB-08-010220-2	57	1.11	9.79	8.08	80.3	
L1010-01	OK-01-010220	58	1.18	9.93	9.31	92.9	
L1010-03	OK-01-010220	59	1.15	9.77	9.4	95.7	
L1010-04	OK-01-010220-1	60	1.19	9.74	9.29	94.7	
L1010-05	OK-01-010220-2	61	1.14	9.84	9.22	92.9	
L1011-01	MC-010220-0-1-GRAB-B3	40	1.13	9.62	8.22	83.5	
L1011-02	MC-010220-0-2-COMP-B3	41	1.16	9.87	8.02	78.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 1/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 01/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107005

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1011-03	MC-010220-4-5-GRAB-B3	42	1.17	9.62	8.48	86.5	
L1011-04	MC-010220-2-5-COMP-B3	43	1.11	9.52	8.07	82.8	
L1011-05	MC-010220-0-1-GRAB-B7	44	1.15	9.94	8.71	86.0	
L1011-06	MC-010220-0-2-COMP-B7	45	1.16	9.64	8.22	83.3	
L1011-07	MC-010220-4-5-GRAB-B7	46	1.17	9.92	8.63	85.3	
L1011-08	MC-010220-2-5-COMP-B7	47	1.12	9.55	7.22	72.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$