



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 1/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:00
Out Date: 01/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107077

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1029-03	X068-01-I	33	1.15	9.64	8.13	82.2	
L1029-06	X068-02-I	34	1.16	9.96	8.2	80.0	
L1029-09	X068-03-I	35	1.18	9.9	7.66	74.3	
L1029-13	X068-04-I	36	1.13	9.57	7.4	74.3	
L1029-16	X068-05-I	37	1.18	9.9	7.24	69.5	
L1046-01	C1Y-SB-04-0-2	1	1.12	9.93	8.97	89.1	
L1046-02	C1Y-SB-104-0-2	2	1.13	9.82	9.04	91.0	
L1046-03	C1Y-SB-04-3-4	3	1.14	9.59	8.23	83.9	
L1046-04	C1Y-SB-05-0-2	4	1.15	9.94	9.18	91.4	
L1046-05	C1Y-SB-05-3-4	5	1.12	9.61	8.81	90.6	
L1046-06	C1Y-SB-06-0-2	6	1.18	9.65	8.42	85.5	
L1046-07	C1Y-SB-06-3-4	7	1.16	9.97	8.42	82.4	
L1046-08	C1Y-SB-07-0-2	8	1.12	9.94	8.5	83.7	
L1046-09	C1Y-SB-07-3-4	9	1.14	9.94	8.58	84.5	
L1046-12	L1046-01MS	10	1.12	9.93	8.97	89.1	
L1046-13	L1046-01MSD	11	1.12	9.93	8.97	89.1	
L1051-01	SB-4G- (16)	12	1.14	9.93	8.29	81.3	
L1051-02	SB-4G- (16-18)	13	1.13	9.81	8.47	84.6	
L1051-03	SB-4G- (18-20)	14	1.14	9.58	8.02	81.5	
L1051-04	SB-4H- (14.5-16.5)	15	1.15	9.92	8.64	85.4	
L1051-05	SB-4F- (22-24)	16	1.14	9.55	8.03	81.9	
L1051-06	SB-3F- (8.5-10.5)	17	1.15	9.96	8.31	81.3	
L1051-07	SB-3G- (11-13)	18	1.14	9.72	7.88	78.6	
L1051-08	SB-3I- (8.5-10.5)	19	1.15	9.95	8.15	79.5	
L1051-09	SB-3H- (8.5-10.5)	20	1.16	9.79	7.8	76.9	
L1051-10	L1051-09MS	21	1.16	9.79	7.8	76.9	
L1051-11	L1051-09MSD	22	1.16	9.79	7.8	76.9	
L1054-02	C1Y-SB-02-0-2	48	1.12	9.53	8.29	85.3	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 1/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:00
Out Date: 01/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107077

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1054-03	C1Y-SB-02-3-4	49	1.19	9.66	8.56	87.0	
L1054-04	C1Y-SB-03-0-2	50	1.19	9.8	8.89	89.4	
L1054-05	C1Y-SB-03-3-4	51	1.15	9.95	9.39	93.6	
L1054-06	C1Y-SB-01-0-2	52	1.1	9.71	8.43	85.1	
L1054-07	C1Y-SB-01-3-4	53	1.18	9.64	8.24	83.5	
L1056-01	WP022SB446_200106_35-36	23	1.15	9.61	8.85	91.0	
L1056-02	WP022SB447_200106_33-34	24	1.18	9.97	8.71	85.7	
L1056-03	WP022SB447_200106_35-36	25	1.14	9.6	8.48	86.8	
L1056-04	WP022SB449_200107_28-30	26	1.19	9.77	8.57	86.0	
L1056-05	WP022SB449_200107_32-34	27	1.16	9.67	8.3	83.9	
L1056-06	WP022SB449_200107_36-38	28	1.17	9.98	8.69	85.4	
L1056-07	WP022SB449_200107_36-38_FD	29	1.19	9.75	7.39	72.4	
L1056-08	WP022SB448_200107_36-38	30	1.13	9.51	8.27	85.2	
L1056-09	WP022SB450_200107_30-32	31	1.15	9.82	8.53	85.1	
L1056-10	WP022SB450_200107_35-36	32	1.1	9.6	8.39	85.8	
L1060-01	TP-1	38	1.14	9.94	8.18	80.0	
L1060-02	TP-1A	39	1.18	9.91	8.27	81.2	
L1060-03	TP-1B	40	1.18	9.91	8.27	81.2	
L1060-04	TP-1	41	1.00	2.00	2.00	100.0	sub
L1060-07	TP-2	42	1.19	9.98	8.39	81.9	
L1060-08	TP-2A	43	1.16	9.62	7.95	80.3	
L1060-09	TP-2B	44	1.16	9.62	7.95	80.3	
L1060-13	TP-3	45	1.15	9.72	7.98	79.7	
L1060-14	TP-3A	46	1.16	9.79	8.44	84.4	
L1060-15	TP-3B	47	1.16	9.79	8.44	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$