



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 1/14/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 01/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:57
Out Date: 01/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107141

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1109-01	GASG2	1	1.12	9.62	7.84	79.1	
L1109-02	GASG2MS	2	1.12	9.62	7.84	79.1	
L1109-03	GASG2MSD	3	1.12	9.62	7.84	79.1	
L1109-04	GASG3	4	1.16	9.74	8.44	84.8	
L1109-05	GASG4	5	1.18	9.84	8.8	88.0	
L1109-06	GASG5	6	1.18	9.9	9.57	96.2	
L1109-07	GASG6	7	1.15	9.57	9.25	96.2	
L1109-08	GASG7	8	1.13	9.6	9.27	96.1	
L1109-09	GASG8	9	1.14	9.97	8.5	83.4	
L1109-10	GASG9	10	1.15	9.81	8.77	88.0	
L1109-11	GASH0	11	1.15	9.96	9.01	89.2	
L1109-12	GASH1	12	1.13	9.89	8.97	89.5	
L1109-13	GASH2	13	1.14	9.95	9.64	96.5	
L1109-14	GASH3	14	1.13	9.92	8.34	82.0	
L1109-15	GASH4	15	1.19	9.97	9.18	91.0	
L1109-16	GASH5	16	1.14	9.79	8.39	83.8	
L1109-17	GASH6	17	1.13	9.58	8.32	85.1	
L1109-18	GASH7	18	1.14	9.57	9.31	96.9	
L1109-19	GASH8	19	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1109-20	GASH9	20	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1109-26	VHBLK02	21	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$