



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/21/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 01/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107260

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1004-13	BU-03-012020	39	1.14	9.84	8.87	88.9	
L1004-14	BU-03-012020	40	1.15	9.95	9.16	91.0	
L1189-01	SB-1-13.5-14.0	1	1.13	9.68	8.87	90.5	
L1189-02	SB-2-13.5-14.0	2	1.13	9.87	9.44	95.1	
L1189-03	SB-3-13.5-14.0	3	1.14	9.58	8.78	90.5	
L1189-04	SB-4-13.5-14.0	4	1.15	9.55	8.34	85.6	
L1189-05	SB-5-13.5-14.0	5	1.13	9.59	8.82	90.9	
L1189-06	SB-6-13.5-14.0	6	1.15	9.63	8.94	91.9	
L1189-07	SB-7-13.5-14.0	7	1.14	9.9	8.69	86.2	
L1189-08	SB-8-13.5-14.0	8	1.15	9.93	8.6	84.9	
L1190-01	TSP-02	9	1.15	9.52	7.64	77.5	
L1190-02	TSP-02	10	1.15	9.79	6.77	65.0	
L1196-01	400114005	11	1.15	9.53	1.67	6.2	SLUDGE SAMPLE
L1197-01	400115768	12	1.13	9.57	3.04	22.6	SLUDGE SAMPLE
L1198-01	400114060	13	1.14	9.98	1.56	4.8	SLUDGE SAMPLE
L1199-01	400115003	14	1.14	9.96	1.74	6.8	SLUDGE SAMPLE
L1200-01	400116710	15	1.14	9.6	2.58	17.0	SLUDGE SAMPLE
L1201-01	400116126	16	1.13	9.5	1.4	3.2	SLUDGE SAMPLE
L1205-01	EP-6	17	1.14	9.75	9.57	97.9	
L1205-02	EP-6	18	1.14	9.89	9.71	97.9	
L1205-03	EP-7	19	1.15	9.82	9.4	95.2	
L1205-04	EP-7	20	1.14	9.83	9.47	95.9	
L1205-05	EP-8	21	1.13	9.76	9.38	95.6	
L1205-06	EP-8	22	1.14	9.54	9.25	96.5	
L1205-07	EP-1	23	1.13	9.75	9.26	94.3	
L1205-08	EP-1	24	1.13	9.6	9.05	93.5	
L1205-09	EP-2	25	1.15	9.63	8.85	90.8	
L1205-10	EP-2	26	1.14	9.97	9.53	95.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/21/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 01/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107260

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1205-11	EP-3	27	1.15	9.69	9.48	97.5	
L1205-12	EP-3	28	1.14	9.92	9.54	95.7	
L1205-13	EP-4	29	1.15	9.62	9.35	96.8	
L1205-14	EP-4	30	1.15	9.55	9.32	97.3	
L1205-15	EP-5	31	1.15	9.71	9.08	92.6	
L1205-16	EP-5	32	1.14	9.62	9.02	92.9	
L1210-01	Q012-SB1	33	1.16	6.34	5.75	88.6	
L1210-02	Q012-SB2	34	1.16	5.16	4.94	94.5	
L1210-03	Q012-SB3	35	1.16	2.78	2.7	95.1	
L1210-04	Q012-SB4	36	1.13	7.04	6.12	84.4	
L1210-05	Q012-SB5	37	1.13	6.94	6.1	85.5	
L1210-06	Q012-WC1	38	1.14	9.83	9.06	91.1	
L1211-01	B1-0-15	41	1.18	9.98	9.27	91.9	
L1211-02	B1-0-15	42	1.18	9.98	9.27	91.9	
L1211-03	B1-15-30	43	1.14	9.94	9.32	93.0	
L1211-04	B1-15-30	44	1.14	9.94	9.32	93.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$