



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/22/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:24
In Date: 01/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 06:45
Out Date: 01/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107284

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1213-01	TP-1	14	1.15	9.65	8.61	87.8	
L1213-03	TP-2	15	1.14	9.73	8.5	85.7	
L1214-01	K500-1	1	1.13	9.54	8.93	92.7	
L1214-02	K500-2	2	1.15	9.6	9.02	93.1	
L1214-03	K500-3	3	1.13	9.72	9.02	91.9	
L1214-04	K500-4	4	1.13	9.74	8.98	91.2	
L1214-05	K500-5	5	1.15	9.91	9.22	92.1	
L1215-01	RBR-200055	16	1.14	9.58	8.28	84.6	
L1215-02	RBR-200031	17	1.15	9.92	8.35	82.1	
L1216-01	SB-01-6.5-7.0	6	1.15	9.85	8.93	89.4	
L1216-02	SB-01-COMP	7	1.14	9.73	8.53	86.0	
L1216-03	SB-02-2.5-3.0	8	1.13	9.8	8.12	80.6	
L1216-04	SB-02-COMP	9	1.15	9.87	8.15	80.3	
L1216-05	SB-03-14.5-15.0	10	1.15	9.58	9.26	96.2	
L1216-06	SB-03-COMP	11	1.14	9.63	8.98	92.3	
L1216-07	SB-04-9.5-10.0	12	1.14	9.66	9.22	94.8	
L1216-08	SB-04-COMP	13	1.13	9.63	9.08	93.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$