



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 01/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107314

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1007-01	JC-03-012320-A	50	1.12	9.8	8.95	90.2	
L1007-02	JC-03-012320-A	51	1.14	9.96	9.48	94.6	
L1007-03	JC-03-012320-A-1	52	1.14	9.6	9.11	94.2	
L1007-04	JC-03-012320-A-2	53	1.13	9.96	9.32	92.8	
L1007-05	JC-03-012320-A-3	54	1.14	9.52	8.8	91.4	
L1007-06	JC-03-012320-A-4	55	1.17	9.74	8.65	87.3	
L1007-07	JC-03-012320-A-5	56	1.14	9.56	9.03	93.7	
L1007-08	JC-03-012320-A-6	57	1.17	9.75	8.54	85.9	
L1007-09	JC-03-012320-A-7	58	1.17	9.9	8.56	84.7	
L1007-11	JC-03-012320-B	59	1.13	9.51	8.77	91.2	
L1007-12	JC-03-012320-B	60	1.18	9.93	9.3	92.8	
L1007-13	JC-03-012320-1	61	1.13	9.54	8.68	89.8	
L1007-14	JC-03-012320-2	62	1.13	9.77	8.88	89.7	
L1007-15	JC-03-012320-3	63	1.14	9.98	9.61	95.8	
L1007-16	JC-03-012320-4	64	1.12	9.94	9.24	92.1	
L1007-17	JC-03-012320-5	65	1.11	9.61	9.13	94.4	
L1007-19	JC-03-012320-C	66	1.18	9.96	9.46	94.3	
L1007-20	JC-03-012320-C	67	1.18	9.93	9.53	95.4	
L1007-21	JC-03-012320-1	68	1.19	9.77	8.93	90.2	
L1007-22	JC-03-012320-2	69	1.13	9.9	9.28	92.9	
L1007-23	JC-03-012320-3	70	1.1	9.78	8.9	89.9	
L1007-24	JC-03-012320-4	71	1.14	9.66	8.77	89.6	
L1007-25	JC-03-012320-5	72	1.15	9.95	9.29	92.5	
L1016-10	SU-04-012320	73	1.13	9.67	8.91	91.1	
L1016-11	SU-04-012320	74	1.14	9.72	9.17	93.6	
L1016-12	SU-04-012320-1	75	1.1	9.79	9.09	91.9	
L1016-13	SU-04-012320-2	76	1.15	9.98	9.4	93.4	
L1161-04	LOD-MDL-SOIL-04-QT1-2020	1	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 01/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107314

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1161-05	LOQ-SOIL-05-QT1-2020	2	1.00	2.00	2.00	100.0	
L1161-06	MDL-SOIL-06-QT1-2020	3	1.00	2.00	2.00	100.0	
L1234-01	NS-TP	4	1.15	9.64	8.01	80.8	
L1234-02	NS-TP	5	1.13	9.65	8.3	84.2	
L1234-03	NS-TP	6	1.15	9.64	8.01	80.8	
L1236-16	40218-S	7	1.15	9.85	5.98	55.5	
L1237-01	PATERSON-UST	77	1.14	9.94	8.53	84.0	
L1238-01	CLIFTON-POLES	78	1.16	9.98	9.31	92.4	
L1239-01	T-1-PE-1	8	1.12	9.55	8.45	87.0	
L1239-02	T-1-PE-2	9	1.13	9.72	8.23	82.7	
L1239-03	T-1-PE-3	10	1.19	9.92	8.54	84.2	
L1239-04	T-1-PE-4	11	1.16	9.95	8.39	82.3	
L1239-05	T-1-PE-5-1	12	1.18	9.72	8.28	83.1	
L1239-06	T-2-PE-1-1	13	1.18	9.98	9.24	91.6	
L1239-07	T-2-PE-2	14	1.14	9.94	8.5	83.6	
L1239-08	T-2-PE-3	15	1.17	9.97	8.76	86.3	
L1239-09	T-2-PE-4	16	1.13	9.87	8.6	85.5	
L1239-10	T-2-PE-5	17	1.13	9.79	8.8	88.6	
L1239-11	T-2-PE-6	18	1.17	9.83	9.51	96.3	
L1239-12	T-2-PE-7	19	1.13	9.94	8.82	87.3	
L1240-01	SB-05-11.5-12.0	30	1.15	9.92	9.22	92.0	
L1240-02	SB-05-COMP	31	1.14	9.8	8.93	90.0	
L1241-01	EP-9	20	1.14	9.61	9.44	98.0	
L1241-02	EP-10	21	1.18	9.55	9.25	96.4	
L1241-03	EP-11	22	1.17	9.95	8.38	82.1	
L1241-04	EP-12	23	1.19	9.51	9.22	96.5	
L1241-05	L1241-04MS	24	1.19	9.51	9.22	96.5	
L1241-06	L1241-04MSD	25	1.19	9.51	9.22	96.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 01/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107314

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1241-07	EP-13	26	1.12	9.74	9.44	96.5	
L1241-08	EP-14	27	1.19	9.97	8.77	86.3	
L1241-09	EP-15	28	1.15	9.73	9.28	94.8	
L1241-10	EP-16	29	1.11	9.72	9.14	93.3	
L1242-01	20200055-COMPOSITE-1	79	1.12	9.88	4.67	40.5	
L1242-02	20200056-COMPOSITE-2	80	1.15	9.96	4.92	42.8	
L1242-03	20200057-COMPOSITE-3	81	1.11	9.92	4.49	38.4	
L1242-04	20200058-COMPOSITE-4	82	1.13	9.58	4.76	43.0	
L1242-05	20200059-COMPOSITE-5	83	1.17	9.5	4.08	34.9	
L1242-06	20200054-FIELD-DUP-COMP	84	1.13	9.96	4.59	39.2	
L1242-07	20200060-COMPOSITE-B6	85	1.13	9.53	4.39	38.8	
L1242-08	20200061-COMPOSITE-B7	86	1.13	9.76	4.52	39.3	
L1242-09	20200062-COMPOSITE-B8	87	1.14	9.75	4.36	37.4	
L1242-10	20200063-AMENDED-COMP-1	88	1.15	9.79	5.27	47.7	
L1242-11	20200064-AMENDED-COMP-2	89	1.15	9.95	5.16	45.6	
L1242-12	20200065-AMENDED-COMP-3	90	1.12	9.63	4.85	43.8	
L1242-13	20200066-AMENDED-COMP-4	91	1.13	9.6	4.81	43.4	
L1242-14	20200067-AMENDED-COMP-5	92	1.17	9.94	5.05	44.2	
L1243-01	CTY-SB-06-1.5-2.5	32	1.16	9.75	8.68	87.5	
L1243-02	CTY-SB-01-1-2	33	1.13	9.76	8.66	87.3	
L1243-03	CTY-SB-101-1-2	34	1.14	9.53	8.54	88.2	
L1243-04	CTY-SB-04-1-2	35	1.13	9.8	8.89	89.5	
L1243-05	CTY-SB-03-1-2	36	1.16	9.93	8.77	86.8	
L1243-06	CTY-SB-02-1-2	37	1.16	9.88	8.53	84.5	
L1243-07	CTY-SB-04-0-5	38	1.17	9.5	8.39	86.7	
L1243-08	CTY-SB-BC-0-3.5	39	1.16	9.83	8.47	84.3	
L1243-09	CTY-SB-01-6-7	40	1.13	9.59	8.5	87.1	
L1243-10	CTY-SB-03-5-6	41	1.14	9.8	8.69	87.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 01/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107314

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1244-01	CTY-WC-S-04	42	1.14	9.61	8.46	86.4	
L1244-02	CTY-WC-S-BC	43	1.15	9.84	8.52	84.8	
L1246-01	SP-1	44	1.13	9.78	8.87	89.5	
L1246-02	SP-1	45	1.14	9.76	8.93	90.4	
L1246-03	SP-2	46	1.14	9.52	8.49	87.7	
L1246-04	SP-2	47	1.16	9.62	8.55	87.4	
L1246-05	SP-3	48	1.15	9.7	8.1	81.3	
L1246-06	SP-3	49	1.13	9.76	8.19	81.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$