



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/29/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 01/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 111
Time OUT: 07:51
Out Date: 01/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107395

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1013-07	PL-05-012820	1	1.13	9.96	9.00	89.1	
L1013-08	PL-05-012820	2	1.12	9.71	8.52	86.1	
L1013-09	PL-05-012820-1	3	1.12	9.78	9.05	91.6	
L1013-10	PL-05-012820-2	4	1.14	9.61	8.75	89.8	
L1013-11	PL-05-012820-3	5	1.15	9.9	8.74	86.7	
L1268-06	17562	6	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L1274-01	BORING-B1-SAMPLE-S1	7	1.13	9.56	9.14	95.0	
L1274-02	BORING-B2-SAMPLE-S1	8	1.14	9.86	9.44	95.2	
L1278-01	TONNELE-AVE-01	9	1.15	9.57	8.61	88.6	
L1287-01	AOC-3-EXC-PIT- (0-5) -A	10	1.13	9.97	8.48	83.1	
L1287-02	AOC-3-EXC-PIT- (5) -A	11	1.14	9.9	8.49	83.9	
L1287-03	AOC-3-EXC-PIT- (0-5) -B	12	1.15	9.6	8.37	85.4	
L1287-04	AOC-3-EXC-PIT- (5) -B	13	1.15	9.54	7.78	79.0	
L1289-01	TP-26-27	14	1.15	7.65	6.37	80.3	
L1289-02	TP-26-27-1.0-1.5	15	1.15	9.93	8.01	78.1	
L1289-03	TP-26-27-2.5-3.0	16	1.15	9.76	7.7	76.1	
L1289-04	TP-26-27-4.0-4.5	17	1.12	9.94	7.74	75.1	
L1289-05	TP-26-27-COMP	18	1.13	9.56	7.82	79.4	
L1289-06	TP-26-27-COMP	19	1.13	9.56	7.82	79.4	
L1289-07	TP-26-27-COMP	20	1.13	9.56	7.82	79.4	
L1289-09	TP-29-30	21	1.13	5.96	5.73	95.2	
L1289-10	TP-29-30-2.0-2.5	22	1.14	9.63	8.77	89.9	
L1289-11	TP-29-30-7.0-7.5	23	1.14	9.98	8.93	88.1	
L1289-12	TP-29-30-11.5-12.0	24	1.14	9.76	8.44	84.7	
L1289-13	TP-29-30-COMP	25	1.13	9.74	8.74	88.4	
L1289-14	TP-29-30-COMP	26	1.13	9.74	8.74	88.4	
L1289-15	TP-29-30-COMP	27	1.13	9.74	8.74	88.4	
L1289-17	TP-32-33	28	1.12	5.67	4.99	85.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/29/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 01/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 111
Time OUT: 07:51
Out Date: 01/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107395

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1289-18	TP-32-33-2.0-2.5	29	1.13	9.54	8.29	85.1	
L1289-19	TP-32-33-7.0-7.5	30	1.13	9.89	8.91	88.8	
L1289-20	TP-32-33-11.5-12.0	31	1.13	9.61	8.23	83.7	
L1289-21	TP-32-33-COMP	32	1.14	9.7	8.56	86.7	
L1289-22	TP-32-33-COMP	33	1.14	9.7	8.56	86.7	
L1289-23	TP-32-33-COMP	34	1.14	9.7	8.56	86.7	
L1289-25	TP-35-34	35	1.14	6.35	4.92	72.6	
L1289-26	TP-35-34-2.0-2.5	36	1.15	9.98	7.68	74.0	
L1289-27	TP-35-34-7.0-7.5	37	1.15	9.62	7.69	77.2	
L1289-28	TP-35-34-11.5-12.0	38	1.14	9.51	7.91	80.9	
L1289-29	TP-35-34-COMP	39	1.12	9.81	8.33	83.0	
L1289-30	TP-35-34-COMP	40	1.12	9.81	8.33	83.0	
L1289-31	TP-35-34-COMP	41	1.12	9.81	8.33	83.0	
L1289-33	TP-3	42	1.15	5.9	3.94	58.7	
L1289-34	TP-3-2.0-2.5	43	1.14	9.92	5.38	48.3	
L1289-35	TP-3-7.0-7.5	44	1.13	9.86	8.53	84.8	
L1289-36	TP-3-11.5-12.0	45	1.12	9.76	7.93	78.8	
L1289-37	TP-3-COMP	46	1.13	9.98	8.29	80.9	
L1289-38	TP-3-COMP	47	1.13	9.98	8.29	80.9	
L1289-39	TP-3-COMP	48	1.13	9.98	8.29	80.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$