



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/31/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:25
In Date: 01/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 01/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107442

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1324-01	GAT58	1	1.15	9.73	8.04	80.3	
L1324-02	GAT59	2	1.18	9.69	8.22	82.7	
L1324-03	GAT60	3	1.15	9.95	8.74	86.3	
L1324-04	GAT61	4	1.11	9.53	8.4	86.6	
L1324-05	GAT62	5	1.18	9.94	9.62	96.3	
L1324-06	GAT63	6	1.17	9.98	9.24	91.6	
L1324-07	GAT64	7	1.18	9.83	7.99	78.7	
L1324-08	GAT65	8	1.18	9.96	8.62	84.7	
L1324-09	GAT66	9	1.18	9.59	8.56	87.8	
L1324-10	GAT66MS	10	1.18	9.59	8.56	87.8	
L1324-11	GAT66MSD	11	1.18	9.59	8.56	87.8	
L1324-12	GAT67	12	1.15	9.81	8.5	84.9	
L1324-13	GAT68	13	1.18	9.96	9.26	92.0	
L1324-14	GAT69	14	1.18	9.87	8.22	81.0	
L1324-15	GAT70	15	1.17	9.97	9.48	94.4	
L1324-16	GAT71	16	1.18	9.83	8.59	85.7	
L1324-17	GAT72	17	1.17	9.97	9.36	93.1	
L1324-18	GAT73	18	1.13	9.96	9.4	93.7	
L1324-19	GAT74	19	1.13	9.77	9.45	96.3	
L1324-20	VHBLK01	20	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$