



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/4/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 02/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107514

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1308-04	STOCK-PILE-B	1	1.15	9.93	8.26	81.0	
L1347-01	GRAB-01	2	1.13	9.58	8.52	87.5	
L1347-02	COMP-01	3	1.16	9.97	8.8	86.7	
L1360-01	TP-1-TANK-BOTTOM	4	1.13	9.86	9.59	96.9	
L1360-02	TP-2	5	1.13	9.83	9.55	96.8	
L1360-03	TP-3	6	1.12	9.8	9.52	96.8	
L1360-04	TP-4	7	1.13	9.81	9.49	96.3	
L1360-06	TP-5	8	1.13	9.91	9.52	95.6	
L1363-01	HR-06-020320	9	1.13	9.89	8.87	88.4	
L1363-02	HR-06-020320	10	1.14	9.59	8.51	87.2	
L1363-03	HR-06-020320-1	11	1.15	9.63	8.7	89.0	
L1363-04	HR-06-020320-2	12	1.14	9.94	8.96	88.9	
L1368-01	EO-01-020320	13	1.14	9.96	9.02	89.3	
L1368-02	EO-01-020320	14	1.14	9.93	8.85	87.7	
L1368-03	EO-01-020320-1	15	1.15	9.85	8.93	89.4	
L1368-04	EO-01-020320-2	16	1.13	9.9	8.61	85.3	
L1368-05	EO-01-020320-3	17	1.15	9.97	9.09	90.0	
L1376-01	02-03-A	18	1.00	2.00	2.00	100.0	100% SOLIDS
L1376-02	02-03-B	19	1.00	2.00	2.00	100.0	100% SOLIDS
L1376-03	02-03-C	20	1.00	2.00	2.00	100.0	100% SOLIDS
L1376-04	OILY-DEBRIS	21	1.00	2.00	2.00	100.0	100% SOLIDS
L1379-01	M030-1A	22	1.15	9.95	8.78	86.7	
L1379-02	M030-1B	23	1.15	9.9	8.65	85.7	
L1379-03	M030-1C	24	1.14	9.98	8.02	77.8	
L1379-04	M030-2A	25	1.15	9.92	8.48	83.6	
L1379-05	M030-2B	26	1.14	9.75	8.56	86.2	
L1379-06	M030-2C	27	1.14	9.68	7.81	78.1	
L1379-07	M030-3A	28	1.13	9.87	6.65	63.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/4/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 02/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107514

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1379-08	M030-3B	29	1.13	9.85	8.56	85.2	
L1379-09	M030-3C	30	1.14	9.61	8.06	81.7	
L1379-10	M030-4A	31	1.14	9.98	8.6	84.4	
L1379-11	M030-4B	32	1.13	9.83	8.6	85.9	
L1379-12	M030-4C	33	1.13	9.86	8.36	82.8	
L1379-13	M030-5A	34	1.14	9.75	4.3	36.7	
L1379-14	M030-6A	35	1.14	9.59	4.7	42.1	
L1379-15	M030-7A	36	1.14	9.95	4.08	33.4	
L1379-16	M030-8A	37	1.12	9.66	4.72	42.2	
L1379-17	M030-9A	38	1.14	9.75	4.1	34.4	
L1379-18	M030-10A	39	1.15	9.66	5.00	45.2	
L1379-19	M030-11A	40	1.13	9.58	4.43	39.1	
L1379-20	M030-12A	41	1.18	9.7	4.02	33.3	
L1380-01	M030-13C	42	1.19	9.64	8.23	83.3	
L1380-02	M030-14C	43	1.13	9.89	8.15	80.1	
L1380-03	M030-15C	44	1.14	9.7	7.52	74.5	
L1380-04	M030-16C	45	1.12	9.92	7.44	71.8	
L1380-05	M030-17C	46	1.13	9.99	8.04	78.0	
L1380-06	M030-18C	47	1.15	9.68	7.4	73.3	
L1380-07	M030-19C	48	1.12	9.97	8.22	80.2	
L1380-08	M030-20C	49	1.13	9.96	8.03	78.1	
L1380-09	M030-21C	50	1.13	9.98	7.56	72.7	
L1380-10	M030-DUP-1	51	1.11	9.66	8.39	85.1	
L1380-11	M030-DUP-2	52	1.12	9.56	8.1	82.7	
L1380-12	M030-DUP-3	53	1.12	9.94	4.19	34.8	
L1380-13	M030-DUP-4	54	1.12	9.52	7.33	73.9	
L1383-01	TP-A	55	1.13	9.65	8.44	85.8	
L1383-02	TP-A-UPPER	56	1.13	9.79	8.66	87.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/4/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 02/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107514

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1383-03	TP-A-LOWER	57	1.13	9.62	7.28	72.4	
L1383-04	TP-A-1	58	1.13	9.65	8.44	85.8	
L1383-05	TP-B	59	1.13	9.89	9.1	91.0	
L1383-06	TP-B-UPPER	60	1.14	9.65	8.89	91.1	
L1383-07	TP-B-LOWER	61	1.14	9.63	8.78	90.0	
L1383-08	TP-B-1	62	1.13	9.89	9.1	91.0	
L1383-09	TP-C	63	1.14	9.5	8.53	88.4	
L1383-10	TP-C-UPPER	64	1.14	9.84	8.08	79.8	
L1383-11	TP-C-LOWER	65	1.17	9.93	8.83	87.4	
L1383-12	TP-C-1	66	1.14	9.5	8.53	88.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$