



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 02/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 02/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107576

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1400-01	CB6G7	1	1.13	9.69	9.58	98.7	
L1400-02	CB6G8	2	1.15	9.82	9.66	98.2	
L1400-03	CB6G9	3	1.14	9.55	9.44	98.7	
L1400-04	CB6H0	4	1.12	9.51	9.36	98.2	
L1400-05	CB6H1	5	1.14	9.77	9.62	98.3	
L1400-06	CB6H2	6	1.14	9.75	9.65	98.8	
L1400-07	CB6H3	7	1.12	9.87	9.7	98.1	
L1400-08	CB6H4	8	1.16	9.82	9.65	98.0	
L1400-09	CB6H5	9	1.12	9.66	9.52	98.4	
L1400-10	CB6H6	10	1.16	9.64	9.53	98.7	
L1400-11	CB6H7	11	1.18	9.81	9.67	98.4	
L1400-12	CB6H8	12	1.17	9.51	9.4	98.7	
L1400-13	CB6H9	13	1.12	9.92	9.8	98.6	
L1400-14	CB6J0	14	1.14	9.74	9.61	98.5	
L1400-15	CB6J1	15	1.13	9.57	9.44	98.5	
L1400-16	CB6J2	16	1.13	9.59	9.45	98.3	
L1400-17	CB6J3	17	1.14	9.9	9.73	98.1	
L1400-18	CB6T5	18	1.13	9.7	9.58	98.6	
L1400-19	CB6T6	19	1.12	9.88	9.75	98.5	
L1400-20	CB6X7	20	1.13	9.66	9.65	99.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$