



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:15
In Date: 02/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:27
Out Date: 02/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107577

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1401-01	CB6J4	1	1.14	9.54	9.44	98.8	
L1401-02	CB6J5	2	1.14	9.6	9.48	98.6	
L1401-03	CB6J6	3	1.12	9.82	9.71	98.7	
L1401-04	CB6J7	4	1.13	9.68	9.58	98.8	
L1401-05	CB6J8	5	1.13	9.52	9.43	98.9	
L1401-06	CB6J9	6	1.13	9.7	9.61	98.9	
L1401-07	CB6K0	7	1.13	9.57	9.52	99.4	
L1401-08	CB6K1	8	1.12	9.64	9.5	98.4	
L1401-09	CB6K2	9	1.13	9.86	9.8	99.3	
L1401-10	CB6K3	10	1.13	9.6	9.5	98.8	
L1401-11	CB6K4	11	1.14	9.92	9.81	98.7	
L1401-12	CB6K5	12	1.14	9.67	9.55	98.6	
L1401-13	CB6K6	13	1.13	9.63	9.51	98.6	
L1401-14	CB6K7	14	1.12	9.7	9.63	99.2	
L1401-15	CB6K8	15	1.12	9.73	9.66	99.2	
L1401-16	CB6T7	16	1.13	9.61	9.53	99.1	
L1401-17	CB6T8	17	1.14	9.6	9.54	99.3	
L1401-18	CB6T9	18	1.15	9.75	9.63	98.6	
L1401-19	CB6W0	19	1.12	9.82	9.72	98.9	
L1401-20	CB6X8	20	1.12	9.63	9.61	99.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$