



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:30
In Date: 02/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 02/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107583

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1402-01	CB6K9	1	1.12	9.67	9.63	99.5	
L1402-02	CB6L0	2	1.12	9.86	9.8	99.3	
L1402-03	CB6L1	3	1.12	9.63	9.53	98.8	
L1402-04	CB6L2	4	1.13	9.69	9.6	98.9	
L1402-05	CB6L3	5	1.14	9.9	9.79	98.7	
L1402-06	CB6L4	6	1.13	9.56	9.44	98.6	
L1402-07	CB6L5	7	1.14	9.62	9.55	99.2	
L1402-08	CB6L6	8	1.13	9.65	9.56	98.9	
L1402-09	CB6L7	9	1.13	9.76	9.63	98.5	
L1402-10	CB6L8	10	1.14	9.73	9.7	99.7	
L1402-11	CB6L9	11	1.13	9.95	9.9	99.4	
L1402-12	CB6M0	12	1.12	9.57	9.51	99.3	
L1402-13	CB6M1	13	1.12	9.7	9.63	99.2	
L1402-14	CB6M2	14	1.15	9.88	9.81	99.2	
L1402-15	CB6M3	15	1.13	9.84	9.78	99.3	
L1402-16	CB6M4	16	1.12	9.51	9.42	98.9	
L1402-17	CB6W1	17	1.14	9.72	9.63	99.0	
L1402-18	CB6W2	18	1.14	9.7	9.6	98.8	
L1402-19	CB6W3	19	1.13	9.65	9.56	98.9	
L1402-20	CB6X9	20	1.15	9.57	9.55	99.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$