



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/14/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 16:30
In Date: 02/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107752

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1514-01	S-1	1	1.14	9.98	8.9	87.8	
L1514-02	S-2	2	1.15	9.95	9.31	92.7	
L1515-01	PE-1	3	1.12	9.57	8.8	90.9	
L1515-02	PE-2	4	1.16	9.82	8.82	88.5	
L1518-01	OILY-STONES	5	1.14	9.94	9.36	93.4	
L1523-01	COMP-1	6	1.14	9.6	8.15	82.9	
L1523-02	COMP-2	7	1.18	9.55	8.23	84.2	
L1523-03	COMP-3	8	1.17	9.83	8.56	85.3	
L1525-01	SB-1A- (7.5-9.5)	9	1.15	9.81	8.09	80.1	
L1525-02	SB-2B- (18.5-20.5)	10	1.14	9.96	8.89	87.9	
L1528-01	FK-GF-02-048-A	11	1.15	9.73	8.91	90.4	
L1528-02	FK-GF-02-048-A	12	1.19	9.56	8.83	91.3	
L1528-03	FK-GF-02-048-B	13	1.19	9.76	8.96	90.7	
L1528-04	FK-GF-02-048-B	14	1.17	9.66	8.82	90.1	
L1529-01	FK-GF-02-049-A	15	1.19	9.71	8.91	90.6	
L1529-02	FK-GF-02-049-A	16	1.15	9.83	9.02	90.7	
L1529-03	FK-GF-02-049-B	17	1.13	9.6	8.81	90.7	
L1529-04	FK-GF-02-049-B	18	1.18	9.82	8.99	90.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$