



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 02/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 02/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107778

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1365-09	CL-02-021420	1	1.12	9.97	9.06	89.7	
L1365-10	CL-02-021420	2	1.13	9.7	8.58	86.9	
L1365-11	CL-02-021420-1	3	1.14	9.84	8.38	83.2	
L1365-12	CL-02-021420-2	4	1.14	9.58	8.22	83.9	
L1365-13	CL-02-021420-3	5	1.13	9.52	8.51	88.0	
L1365-14	CL-02-021420-4	6	1.16	9.68	8.51	86.3	
L1365-15	CL-02-021420-5	7	1.16	9.53	9.18	95.8	
L1530-06	IDW-MCM-200213	8	1.16	9.97	7.85	75.9	
L1530-07	IDW-LGAC-200213	9	1.16	9.73	5.5	50.6	
L1530-08	IDW-GS-200213	10	1.16	9.68	7.45	73.8	
L1530-09	IDW-VCAG-200213	11	1.15	9.64	9.28	95.8	
L1531-01	WFD-SB-03-1.5-2.5	12	1.16	9.6	8.57	87.8	
L1531-02	WFD-SB-103-1.5-2.5	13	1.14	9.95	9.45	94.3	
L1531-03	WFD-SB-02-2-3	14	1.15	9.92	8.89	88.3	
L1533-05	SVOC-GPC-BLANK	15	1.00	2.00	2.00	100.0	
L1533-06	PEST-GPC-BLANK	16	1.00	2.00	2.00	100.0	
L1533-07	PEST-GPC-BLANK-SPIKE	17	1.00	2.00	2.00	100.0	
L1533-08	PCB-GPC-BLANK	18	1.00	2.00	2.00	100.0	
L1533-09	PCB-GPC-BLANK-SPIKE	19	1.00	2.00	2.00	100.0	
L1533-10	SVOC-GPC2-BLANK	20	1.00	2.00	2.00	100.0	
L1533-11	PEST-GPC2-BLANK	21	1.00	2.00	2.00	100.0	
L1533-12	PEST-GPC2-BLANK-SPIKE	22	1.00	2.00	2.00	100.0	
L1533-13	PCB-GPC2-BLANK	23	1.00	2.00	2.00	100.0	
L1533-14	PCB-GPC2-BLANK-SPIKE	24	1.00	2.00	2.00	100.0	
L1544-01	2-5-1	25	1.14	9.81	8.64	86.5	
L1544-02	2-5-1	26	1.14	9.81	8.64	86.5	
L1544-03	2-5-1	27	1.14	9.81	8.64	86.5	
L1544-05	2-4-1	28	1.13	9.9	8.72	86.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 02/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 02/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107778

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1544-06	2-4-1	29	1.13	9.9	8.72	86.5	
L1544-07	2-4-1	30	1.13	9.9	8.72	86.5	
L1545-01	43ST-B	65	1.15	9.55	8.71	90.0	
L1545-02	43ST-A	66	1.12	9.72	8.76	88.8	
L1546-01	INDIAN-ISLAND-(1)	31	1.16	9.83	8.94	89.7	
L1546-02	INDIAN-ISLAND-(1)	32	1.14	9.98	9.14	90.5	
L1546-03	INDIAN-ISLAND-(2)	33	1.16	9.89	8.97	89.5	
L1546-04	INDIAN-ISLAND-(2)	34	1.16	9.66	8.6	87.5	
L1546-05	INDIAN-ISLAND-(3)	35	1.16	9.51	8.24	84.8	
L1546-06	INDIAN-ISLAND-(3)	36	1.14	9.67	8.93	91.3	
L1546-07	INDIAN-ISLAND-(4)	37	1.16	9.71	8.15	81.8	
L1546-08	INDIAN-ISLAND-(4)	38	1.15	9.74	8.22	82.3	
L1546-09	INDIAN-ISLAND-(5)	39	1.13	9.87	8.6	85.5	
L1546-10	INDIAN-ISLAND-(5)	40	1.13	9.96	8.99	89.0	
L1546-11	INDIAN-ISLAND-(6)	41	1.16	9.56	8.2	83.8	
L1546-12	INDIAN-ISLAND-(6)	42	1.16	9.65	8.36	84.8	
L1546-13	INDIAN-ISLAND-(7)	43	1.16	9.77	8.25	82.3	
L1546-14	INDIAN-ISLAND-(7)	44	1.16	9.98	8.59	84.2	
L1546-15	INDIAN-ISLAND-(8)	45	1.15	9.68	8.26	83.4	
L1546-16	INDIAN-ISLAND-(8)	46	1.16	9.85	8.42	83.5	
L1546-17	INDIAN-ISLAND-(9)	47	1.13	9.97	8.59	84.4	
L1546-18	INDIAN-ISLAND-(9)	48	1.16	9.81	8.55	85.4	
L1546-19	INDIAN-ISLAND-(10)	49	1.16	9.88	8.34	82.3	
L1546-20	INDIAN-ISLAND-(10)	50	1.16	9.9	8.5	84.0	
L1546-21	INDIAN-ISLAND-(11)	51	1.15	9.96	8.69	85.6	
L1546-22	INDIAN-ISLAND-(11)	52	1.16	9.79	8.63	86.6	
L1546-23	INDIAN-ISLAND-(12)	53	1.15	9.94	8.13	79.4	
L1546-24	INDIAN-ISLAND-(12)	54	1.16	9.98	8.07	78.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 02/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 02/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107778

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1547-01	MARCY-GREEN-BH-01-A	67	1.18	9.93	8.09	79.0	
L1547-02	MARCY-GREEN-BH-01-B	68	1.19	9.77	8.74	88.0	
L1547-03	MARCY-GREEN-BH-01-C	69	1.17	9.57	8.82	91.1	
L1547-04	MARCY-GREEN-BH-01-A	70	1.11	9.6	7.29	72.8	
L1547-05	MARCY-GREEN-BH-01-B	71	1.17	9.97	9.00	89.0	
L1547-06	MARCY-GREEN-BH-01-C	72	1.13	9.91	9.08	90.5	
L1547-07	MARCY-GREEN-BH-02-A	73	1.15	9.96	9.22	91.6	
L1547-08	MARCY-GREEN-BH-02-B	74	1.11	9.75	8.92	90.4	
L1547-09	MARCY-GREEN-BH-02-C	75	1.16	9.97	9.00	89.0	
L1547-10	MARCY-GREEN-BH-02-A	76	1.13	9.89	8.87	88.4	
L1547-11	MARCY-GREEN-BH-02-B	77	1.13	9.72	8.8	89.3	
L1547-12	MARCY-GREEN-BH-02-C	78	1.16	9.92	8.77	86.9	
L1547-13	MARCY-GREEN-BH-03-A	79	1.14	9.66	8.12	81.9	
L1547-14	MARCY-GREEN-BH-03-B	80	1.17	9.85	8.94	89.5	
L1547-15	MARCY-GREEN-BH-03-C	81	1.13	9.88	8.66	86.1	
L1547-16	MARCY-GREEN-BH-03-A	82	1.14	9.96	8.34	81.6	
L1547-17	MARCY-GREEN-BH-03-B	83	1.19	9.72	8.24	82.6	
L1547-18	MARCY-GREEN-BH-03-C	84	1.19	9.71	8.85	89.9	
L1547-19	MARCY-GREEN-BH-04-A	85	1.17	9.6	9.11	94.2	
L1547-20	MARCY-GREEN-BH-04-B	86	1.19	9.98	9.57	95.3	
L1547-21	MARCY-GREEN-BH-04-C	87	1.15	9.96	9.36	93.2	
L1547-22	MARCY-GREEN-BH-04-A	88	1.14	9.97	9.46	94.2	
L1547-23	MARCY-GREEN-BH-04-B	89	1.1	9.98	9.05	89.5	
L1547-24	MARCY-GREEN-BH-04-C	90	1.17	9.96	9.15	90.8	
L1548-01	FK-PS-01-063	55	1.15	9.84	8.56	85.3	
L1548-02	FK-PS-01-063	56	1.16	9.58	8.36	85.5	
L1551-01	FK-SF-01-007-A	57	1.16	9.61	8.44	86.2	
L1551-02	FK-SF-01-007-A	58	1.16	9.54	8.43	86.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 02/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 02/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107778

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1551-03	FK-SF-01-007-B	59	1.16	9.59	8.38	85.6	
L1551-04	FK-SF-01-007-B	60	1.15	9.68	8.5	86.2	
L1552-01	FK-SF-01-008-A	61	1.16	9.54	8.37	86.0	
L1552-02	FK-SF-01-008-A	62	1.16	9.81	8.7	87.2	
L1552-03	FK-SF-01-008-B	63	1.14	9.66	8.41	85.3	
L1552-04	FK-SF-01-008-B	64	1.16	9.69	8.58	87.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$