



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/20/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 02/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 02/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107846

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1576-01	WFD-SB-010203-0-8	1	1.15	9.95	8.75	86.4	
L1576-02	WFD-SB-1010203-0-8	2	1.19	9.88	8.71	86.5	
L1578-01	WFD-WC-S-SB-010203	3	1.16	9.78	8.68	87.2	
L1581-02	40712	4	1.16	9.77	8.41	84.2	
L1583-02	13192	5	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L1587-01	0170-COMP-SL-A	6	1.12	9.83	4.45	38.2	
L1587-02	0171-COMP-SL-B	7	1.14	9.97	4.96	43.3	
L1587-03	0172-COMP-SL-C	8	1.17	9.95	6.34	58.9	
L1589-01	SB3-9.5-10.0	9	1.12	9.93	8.59	84.8	
L1589-02	SB4-4.5-5.0	10	1.13	9.56	9.36	97.6	
L1589-03	SB5-4.5-5.0	11	1.15	9.98	9.48	94.3	
L1589-04	SB6-3.0-3.5	12	1.17	9.68	9.46	97.4	
L1589-05	SB7-3.0-3.5	13	1.18	9.81	9.48	96.2	
L1589-06	SB8-2.0-2.5	14	1.18	9.64	9.38	96.9	
L1592-01	FK-BPM-04-005-ABCD	15	1.16	9.68	9.22	94.6	
L1592-02	FK-BPM-04-005-ABCD	16	1.17	9.61	8.58	87.8	
L1595-01	FK-BPM-01-C2-C3-C2A-C3A	17	1.16	9.9	9.00	89.7	
L1595-02	FK-BPM-01-C2-C3-C2A-C3A	18	1.19	9.76	8.86	89.5	
L1596-01	FK-PS-01-064	19	1.12	9.57	8.09	82.5	
L1596-02	FK-PS-01-064	20	1.18	9.71	8.49	85.7	
L1599-01	FK-BPM-03-011-ABCD	21	1.11	9.94	8.87	87.9	
L1599-02	FK-BPM-03-011-ABCD	22	1.13	9.77	9.3	94.6	
L1600-01	TP-B-S1	23	1.13	9.81	8.61	86.2	
L1600-02	TP-B-S2	24	1.1	9.84	8.58	85.6	
L1600-03	TP-B-S	25	1.15	9.59	8.34	85.2	
L1600-05	TP-B-D1	26	1.14	9.56	8.76	90.5	
L1600-06	TP-B-D2	27	1.15	9.83	8.06	79.6	
L1600-07	TP-B-D	28	1.18	9.91	8.82	87.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/20/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 02/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 02/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107846

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1601-01	B-2	29	1.15	9.95	8.21	80.2	
L1601-02	B-1	30	1.17	9.81	8.09	80.1	
L1601-03	SW-1	31	1.17	9.98	8.33	81.3	
L1601-04	SW-2	32	1.15	9.93	8.26	81.0	
L1601-05	SW-4	33	1.13	9.86	8.1	79.8	
L1601-06	SW-3	34	1.15	9.78	8.15	81.1	
L1601-07	B-3	35	1.15	9.69	7.98	80.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$