



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/21/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 02/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 02/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107863

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1609-13	MX2L79	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-14	MX2L79D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-15	MX2L79S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-16	MX2L80	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-17	MX2L80D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-18	MX2L80S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-19	MX2L81	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-20	MX2L81D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-21	MX2L81S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-22	MX2L82	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-23	MX2L82D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
L1609-24	MX2L82S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$