



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/28/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 02/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:47
Out Date: 02/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB107974

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1690-01	MGAW29	1	1.14	9.72	9.45	96.9	
L1690-02	MGAW30	2	1.19	9.89	9.43	94.7	
L1690-03	MGAW31	3	1.14	9.83	9.78	99.4	
L1690-04	MGAW32	4	1.12	9.63	9.57	99.3	
L1690-05	MGAW33	5	1.17	9.91	9.82	99.0	
L1690-06	MGAW34	6	1.13	9.51	9.38	98.4	
L1690-07	MGAW35	7	1.13	9.52	9.42	98.8	
L1690-08	MGAW36	8	1.14	9.5	8.88	92.6	
L1690-09	MGAW37	9	1.18	9.79	9.41	95.6	
L1690-10	MGAW38	10	1.12	9.84	9.29	93.7	
L1690-11	MGAW39	11	1.17	9.86	8.63	85.8	
L1690-12	MGAW40	12	1.18	9.92	9.41	94.2	
L1690-13	MGAW40D	13	1.18	9.92	9.41	94.2	
L1690-14	MGAW40S	14	1.18	9.92	9.41	94.2	
L1690-15	MGAW41	15	1.1	9.89	9.38	94.2	
L1690-16	MGAW42	16	1.17	9.68	9.02	92.2	
L1690-17	MGAW43	17	1.17	9.5	8.83	92.0	
L1690-18	MGAW44	18	1.18	9.88	9.2	92.2	
L1690-19	MGAW45	19	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$