



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/28/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 02/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 02/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108000

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1677-01	C0CC5	1	1.14	9.62	6.67	65.2	
L1677-02	C0CC6	2	1.16	9.98	7.08	67.1	
L1677-03	C0CC7	3	1.13	9.92	6.72	63.6	
L1677-04	C0CC8	4	1.16	9.94	7.32	70.2	
L1677-05	C0CC9	5	1.17	9.78	6.29	59.5	
L1677-08	VHBLK02	6	1.00	2.00	2.00	100.0	VHBLK
L1677-09	C0CB9	7	1.17	9.96	7.17	68.3	
L1677-10	C0CC0	8	1.19	9.83	7.36	71.4	
L1677-11	C0CC1	9	1.15	9.96	7.2	68.7	
L1677-12	C0CC2	10	1.1	9.91	6.62	62.7	
L1677-13	C0CD1	11	1.19	9.79	7.88	77.8	
L1677-14	C0CD2	12	1.18	9.57	6.92	68.4	
L1677-15	C0CD3	13	1.17	9.98	6.49	60.4	
L1677-16	C0CD8	14	1.13	9.61	6.09	58.5	
L1677-19	C0CB6	15	1.15	9.91	7.79	75.8	
L1677-20	C0CB7	16	1.14	9.94	6.86	65.0	
L1677-21	C0CB8	17	1.12	9.95	7.46	71.8	
L1677-22	C0CB8MS	18	1.12	9.95	7.46	71.8	
L1677-23	C0CB8MSD	19	1.12	9.95	7.46	71.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$