



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 02/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108018

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1369-01	PL-02-022820	45	1.17	9.89	9.15	91.5	
L1369-02	PL-02-022820	46	1.14	9.98	9.53	94.9	
L1369-03	PL-02-022820-1	47	1.18	9.63	9.2	94.9	
L1369-04	PL-02-022820-2	48	1.18	9.94	9.33	93.0	
L1369-05	PL-02-022820-3	49	1.13	9.9	9.71	97.8	
L1372-06	SU-03-022820	50	1.17	9.96	9.3	92.5	
L1372-07	SU-03-022820	51	1.17	9.63	8.54	87.1	
L1372-08	SU-03-022820-1	52	1.17	9.96	8.79	86.7	
L1372-09	SU-03-022820-2	53	1.1	9.83	9.34	94.4	
L1715-01	S-1-C	1	1.18	9.97	9.22	91.5	
L1715-02	PE-1	2	1.13	9.88	8.47	83.9	
L1715-03	S-2-C	3	1.16	9.96	9.01	89.2	
L1715-04	PE-2	4	1.17	9.65	8.59	87.5	
L1717-05	PLANK-MANHOLE	5	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L1723-01	B3-(1.5-2)	6	1.18	9.98	8.43	82.4	
L1723-02	B2-(1.5-2)	7	1.14	9.95	9.03	89.6	
L1723-03	WC	30	1.12	9.95	7.95	77.3	
L1727-01	PCB-022720-SR-01	8	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L1727-02	PCB-022720-SR-02	9	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L1727-03	PCB-022720-SR-03	10	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L1731-01	OILY-STONE-2-28-2020	54	1.00	2.00	2.00	100.0	OILY-STONE
L1733-05	SVOC-GPC-BLANK	11	1.00	2.00	2.00	100.0	
L1733-06	PEST-GPC-BLANK	12	1.00	2.00	2.00	100.0	
L1733-07	PEST-GPC-BLANK-SPIKE	13	1.00	2.00	2.00	100.0	
L1733-08	PCB-GPC-BLANK	14	1.00	2.00	2.00	100.0	
L1733-09	PCB-GPC-BLANK-SPIKE	15	1.00	2.00	2.00	100.0	
L1733-10	SVOC-GPC2-BLANK	16	1.00	2.00	2.00	100.0	
L1733-11	PEST-GPC2-BLANK	17	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 02/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108018

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1733-12	PEST-GPC2-BLANK-SPIKE	18	1.00	2.00	2.00	100.0	
L1733-13	PCB-GPC2-BLANK	19	1.00	2.00	2.00	100.0	
L1733-14	PCB-GPC2-BLANK-SPIKE	20	1.00	2.00	2.00	100.0	
L1734-01	RW-1-20200228-WCS	21	1.00	2.00	2.00	100.0	WIPE SAMPLE,100%SOLIDS
L1734-02	RW-2-20200228-WCS	22	1.00	2.00	2.00	100.0	WIPE SAMPLE,100%SOLIDS
L1734-03	RW-M-20200228-WCS	23	1.00	2.00	2.00	100.0	WIPE SAMPLE,100%SOLIDS
L1735-01	FK-SF-01-009-A	31	1.14	9.69	9.05	92.5	
L1735-02	FK-SF-01-009-A	32	1.17	9.63	8.93	91.7	
L1735-03	FK-SF-01-009-B	33	1.13	9.5	8.9	92.8	
L1735-04	FK-SF-01-009-B	34	1.18	9.75	9.08	92.2	
L1736-01	FK-SF-01-010-A	35	1.15	9.64	9.04	92.9	
L1736-02	FK-SF-01-010-A	36	1.16	9.54	8.87	92.0	
L1736-03	FK-SF-01-010-B	37	1.17	9.55	8.92	92.5	
L1736-04	FK-SF-01-010-B	38	1.13	9.93	9.26	92.4	
L1737-01	FK-BPM-03-013-ABCD	39	1.18	9.73	8.97	91.1	
L1737-02	FK-BPM-03-013-ABCD	40	1.1	9.81	9.04	91.2	
L1738-01	FK-BPM-01-G0-G1-G0A-G1A	41	1.18	9.69	9.44	97.1	
L1738-02	FK-BPM-01-G0-G1-G0A-G1A	42	1.14	9.6	9.34	96.9	
L1739-01	FK-BPM-04-007-ABCD	43	1.19	9.72	8.82	89.4	
L1739-02	FK-BPM-04-007-ABCD	44	1.16	9.88	9.41	94.6	
L1740-01	GP-06-022720	24	1.13	9.77	7.69	75.9	
L1740-02	GP-07-022720	25	1.13	9.81	8.33	82.9	
L1740-03	GP-08-022720	26	1.16	9.97	9.52	94.9	
L1740-04	GP-09-022720	27	1.17	9.55	7.8	79.1	
L1740-06	GP-12-022720	28	1.17	9.96	9.08	90.0	
L1740-07	GP-13-022720	29	1.12	9.76	8.46	85.0	
L1741-01	TP-1	55	1.16	9.79	8.86	89.2	
L1741-02	TP-1	56	1.16	9.79	8.86	89.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 02/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108018

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1741-03	TP-2	57	1.12	9.54	8.2	84.1	
L1741-04	TP-2	58	1.12	9.54	8.2	84.1	
L1742-01	20C069	59	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100% SOLIDS
L1744-01	TP-3	60	1.15	9.96	9.41	93.8	
L1744-02	TP-2	61	1.18	9.92	9.4	94.1	
L1744-03	TP-1	62	1.1	9.7	8.99	91.7	
L1744-04	TP-4	63	1.13	9.64	9.06	93.2	
L1745-01	MP-2-COMP	64	1.16	9.93	9.01	89.5	
L1745-03	MP-2-2.0-2.5	65	1.19	9.94	9.22	91.8	
L1745-04	MP-2-7.0-7.5	66	1.11	9.77	8.61	86.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$