



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/4/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 17:10
In Date: 03/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 03/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108075

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1721-03	STOCKPILE-SCREENINGS	1	1.18	9.83	6.6	62.7	
L1748-01	AU-06-030320	2	1.18	9.78	8.83	89.0	
L1748-02	AU-06-030320	3	1.19	9.63	8.69	88.9	
L1748-03	AU-06-030320	4	1.17	9.8	8.94	90.0	
L1759-01	OR-02-030320	5	1.12	9.57	8.9	92.1	
L1759-02	OR-02-030320	6	1.18	9.65	8.79	89.8	
L1759-03	OR-02-030320-1	7	1.16	9.86	9.47	95.5	
L1759-04	OR-02-030320-2	8	1.14	9.94	9.46	94.5	
L1771-01	PIER-122-CORE-1-3-2-20	9	1.15	9.72	4.98	44.7	
L1771-02	PIER-122-CORE-2-3-2-20	10	1.15	9.82	5.02	44.6	
L1772-01	EP-9	11	1.17	9.8	9.34	94.7	
L1772-02	EP-7	12	1.16	9.97	9.42	93.8	
L1772-03	EP-8	13	1.18	9.54	8.39	86.2	
L1774-01	VNJ-223	14	1.15	9.57	8.06	82.1	
L1774-02	VNJ-223	15	1.16	9.64	7.91	79.6	
L1774-03	VNJ-236	16	1.14	9.72	8.84	89.7	
L1774-04	VNJ-236	17	1.17	9.9	8.75	86.8	
L1775-01	CR-1	18	1.19	9.96	8.66	85.2	
L1775-02	CR-1-1.0-1.5	19	1.15	9.68	7.58	75.4	
L1775-03	CR-1-2.0-2.5	20	1.13	9.77	7.73	76.4	
L1775-04	CR-1-2.5-3.0	21	1.1	9.52	7.2	72.4	
L1775-07	CR-2	22	1.19	9.8	8.17	81.1	
L1775-08	CR-2-0.5-1.0	23	1.15	9.84	8.09	79.9	
L1775-09	CR-2-2.0-2.5	24	1.15	9.67	8.49	86.2	
L1775-10	CR-2-2.5-3.0	25	1.18	9.71	8.5	85.8	
L1775-13	CR-3	26	1.19	9.8	8.18	81.2	
L1775-14	CR-3-0.5-1.0	27	1.14	9.71	8.03	80.4	
L1775-15	CR-3-1.5-2.0	28	1.15	9.93	8.5	83.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/4/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 17:10
In Date: 03/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 03/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108075

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1775-16	CR-3-2.5-3.0	29	1.18	9.98	8.52	83.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$